

BACTERIAL MIGRATION THROUGH PUNCTURED SURGICAL GLOVES UNDER REAL SURGICAL CONDITIONS

By: Nils Olaf Hübner, Anna Maria Goerdts, Natalie Stanislawski, Ojan Assadian*, Claus-Dieter Heidecke, Axel Kramer, and Lars Ivo Partecke

BACKGROUND

Intact surgical gloves are a barrier against the transmission of blood-borne pathogens. However, recent research indicates that unnoticed micro-perforations in gloves can develop depending on the length of wear and the nature of the surgical procedure performed. This study aimed to validate previous findings, namely that bacterial transmission can occur from the surgical site through the micro-perforations in surgical gloves.



OBJECTIVES AND METHODS


The modified and standardized Gaschen-Bag method was used to identify bacterial migration from the surgical site through the punctured glove. Over a six-month period, double-gloving surgical gloves were used during colorectal surgeries where inner gloves' microbial content were examined. They were subsequently compared with microorganisms detected by an intraoperative swab.

RESULTS

194 gloves were examined, which consisted of 98 outer gloves and 96 inner gloves that were used in 20 consecutive surgeries. Perforations were found in 10 (10.2%) of outer gloves and 1 (1.04%) of inner gloves. The rate of micro-perforation of the outer glove wearing time represents a range of 20-175 minutes with a median wearing time of 100 minutes.

In a total of six cases, bacterial migration was observed through micro-perforations and could be demonstrated microbiologically, as sampled from 20 operations. Bacterial migration through micro-perforation was seen in 5% of the examined outer gloves and 1% of the examined inner gloves. For the outer layer, the total migration was 50% (n=5), representing a minimum wearing time of 62 minutes, and a median wearing time of 71 minutes.

Gloves, perforations/migration detected and performance of the indicator system



CONCLUSION

The study confirms that micro-perforations in surgical gloves allow bacterial migration to occur in real surgical conditions. The detected frequency with which micro-perforations allowed bacterial migration from the patient to the hands of the surgeon was confirmed. This illustrates that the protective barrier function of such compromised gloves is diminished, particularly over extended periods, allowing bacterial migration to occur due to undetected micro-perforation. This information aligns with previous results.



APPLICATION FOR PRACTICE



Change gloves every 90 minutes especially in lengthy surgeries



Double glove with a dark-colored underglove



Follow best clinical practice and assess the risk prior to any surgeries

* Dr. Assadian has performed paid consulting work for Quantum Management & Service, a paid service provider of Ansell and/or its affiliates.

Note: This clinical summary is written by clinicians at Ansell Healthcare Products LLC. Please refer to the actual study for full text information.

Hübner N-O., Gerdt A-M., Stanislawski N., Assadian O., Heidecke C-D., Kramer A., Partecke LI., Bacterial migration through punctured surgical gloves under real surgical conditions BMC Infectious Diseases 2010, 10 :192

➤ For more information or additional clinical resources, please visit: www.ansell.com/AnsellCARES

Ansell, ® and ™ are owned by Ansell Limited or one of its affiliates. © 2023 Ansell Limited. All rights reserved.