# GAMMEX<sup>®</sup> POWDER-FREE WITH AMT INSTRUCTIONS FOR USE

Available in 26 languages  
10/24

Please click on the box below to select your preferred language.

**North America**  
**Manufactured For/Fabricant**  
Ansell Healthcare Products LLC  
Robb Drive, Reno,  
NV 89523, USA 2301

**UK Responsible Person for Medical and  
Importer for PPE Products**  
Nitritex Ltd  
Ground Floor, 15 Kings Court,  
Willie Snaith Road, Newmarket,  
Suffolk, CB8 7SG, United Kingdom

**Asia Pacific**  
Ansell Global Trading Center Malaysia (Sdn Bhd)  
Prima 6, Prima Avenue  
Block 3512, Jalan Teknokrat 6,  
Cyberjaya 63000 Selangor, Malaysia

**EU & UK**  
**Legal Manufacturer**  
Ansell Healthcare Europe N.V.  
Boulevard International 55  
Brussels, B-1070, Belgium

**Australia**  
**Authorised Sponsor ANZ**  
Ansell Limited  
Level 3, 678 Victoria Street Richmond,  
Victoria, 3121

**Manufactured By:**  
Ansell N.P. Sdn. Bhd.  
Lot 80, Air Keroh Industrial Estate  
Melaka 75450 Malaysia

[www.ansell.com](http://www.ansell.com)

Ansell, <sup>®</sup> and <sup>™</sup> are trademarks owned by Ansell Limited or one of its affiliates. © 2024 Ansell Limited. All Rights Reserved.

# CONTENTS

Please click on the content below to navigate to your selected page.

# SYMBOLS & PICTOGRAMS

EN ISO 21420:2020

(01)

EN ISO 374-5:2016

(02)

EN ISO 374-1:2016+A1:  
2018 Type A, B or C

ABCDEFGHIJKLMNOPS

(03)

(04)

(05)

(06)

(07)

(08)

UA.TR.116

(09)

(10)

(11)

(12)

(13)

(14)

(15)

**USE:** This Instruction for Use is to be used in combination with the specific information that is mentioned on each packaging enclosure. These gloves are designed to protect the hands against chemical splashes and comply with the applicable harmonised EN or EN ISO Standards as shown by the pictograms being mentioned on the packaging. The gloves therefore will provide protection against the specific risks as shown by the pictograms which are defined by these harmonised standards. The gloves are in conformity with the European Regulation 2016/425. Please ensure the gloves are used only for the designated purposes, as explained above. The Summary of Safety and Clinical Performance can be found on [ansell.com](http://ansell.com) within the product downloads area. If a serious incident occurs with the glove please report it to Ansell Healthcare Europe NV and the Competent Authority where it was being used. Available in sizes: 5.5/6/6.5/7/7.5/8/8.5/9

**EXPLANATION OF SYMBOLS & PICTOGRAMS:** (01) **EN ISO 21420:2020 / ISO 15223-1:2021** – Please read the Instructions for Use, prior to using the products, or contact Ansell for more information. If a level X is mentioned under any of the pictograms, this means this test is not applicable and glove is not designed and therefore not to be used for this specific hazard. (02) **EN ISO 374-5:2016** – Protection against bacteria, fungi and virus. Not tested against viruses. (03) **EN ISO 374-1:2016+A1:2018 TYPE A, B OR C** – Protection against chemical hazards: Type A = chemical breakthrough time > 30 minutes against at least 6 chemicals as per list defined in EN ISO 374-1:2016 Type B = chemical breakthrough time > 30 minutes against at least three chemicals as per list defined in EN ISO 374-1:2016 Type C = chemical breakthrough time > 10 minutes against at least one test chemical as per list defined in EN ISO 374-1:2016 (no code underneath the pictogram) A = methanol - B = acetone - C = acetonitrile - D = dichloromethane - E = carbon disulfide - F = toluene - G = diethylamine - H = tetrahydrofuran - I = ethyl acetate - J = n-heptane - K = sodium hydroxide, 40% - L = sulphuric acid, 96 % - M = nitric acid, 65% - N = acetic acid, 99% - O = ammonia, 25% - P = hydrogen peroxide, 30% - S = hydrofluoric acid, 40% - T = formaldehyde, 37%. PPE EU -Type examination certificate (Module B) and EU Conformity to type based on quality assurance of the production process (Module D) by Centexbel, Technologiepark 70, B-9052 Zwijnaarde.

**WARNING!** Chemical resistance data provided, has been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. It can be different if it be used in a mixture. The chemical resistance data may not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Degradation results indicate the change in puncture resistance of the gloves after exposure to the challenge chemical. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Chemical permeation data, as tested per EN 16523-1:2015 test method, and degradation data, tested per EN ISO 374-4:2019 test method, are available upon and/or via [ansell.com](http://ansell.com), through the Ansell product page/downloads criteria/CE-certified Chemical Permeation test data. For more detailed information on the product's performance, please consult Ansell. To obtain the EU-Conformity Declaration to PPE Regulation, please use the link as shown hereunder: [www.ansell.com/regulatory](http://www.ansell.com/regulatory) (04) **CE** – Product is compliant and certified to the requirements of the European Regulation on Personal Protective Equipment 2016/425 and European Medical Device Regulation 2017/745. For compliance as a PPE: The CE mark is followed by a four digit code which refers to the identification number of the Notified Body that is in charge of the PPE category III conformity assessment (Module C2 or D), for products to protect against serious risks. (05) **MEDICAL DEVICE (06) UKCA** = Product is compliant to the requirements of the Regulation 2016/425 on personal protective equipment, as amended to apply in GB and UK Medical Devices Regulation 2002 (SI 2002 No. 618 as amended) PPE Type examination certificate (Module B) and Conformity to type based on quality assurance of the production process (Module D) by Centexbel Belgium (I.D. 0493), Technologiepark 70, B-9052 Zwijnaarde, for EU. For UK: Type-examination certificate (Module B) and Conformity to type based on quality assurance of the production process (Module D) for CE marking are used as the basis for applying a UKCA. For more details, please consult UK Conformity Declaration. For compliance as a Medical device product, the CE and UKCA mark is followed by a four digit code which refers to the identification number of the Notified / Approved Body used For EU: I.D. 2797 – BSI Netherlands. For UK: I.D. 0086 – BSI UK. To obtain the MDR EU or UKCA Declaration of Conformity, please visit Ansell Website Product Page/download section: [www.ansell.com/regulatory](http://www.ansell.com/regulatory) (07) **SINGLE USE ONLY (08) A SINGLE STERILE BARRIER SYSTEM, STERILISED USING IRRADIATION (09) UA.TR.116** = Medical device Class III under Technical Regulation on medical devices, approved by Decree of Cabinet of Ministers of Ukraine No. 753 on October 2nd., 2013 under conformity assessment procedure by the Notified Body "UNI-CERT" LLC, with reference (UA.TR.116) (10) **MADE WITH NATURAL RUBBER LATEX (11) DO NOT USE IF PACKAGE IS DAMAGED (12) KEEP DRY (13) KEEP AWAY FROM SUNLIGHT (14) AVOID OZONE (15) EN 455-3:2022 - CONTAINS POTENTIAL TYPE IV CHEMICAL ALLERGENS.** For more detailed information on the product's performance, please consult Ansell.

**PRECAUTIONS FOR USE:** 1. Before usage, inspect the gloves for any defects or imperfections such as holes, pinholes and tears. If the gloves are ripped or punctured during use, dispose of them immediately. If in doubt, do not use the gloves, get a new pair. 2. Do not reverse the gloves. 3. If the gloves are used against chemicals, it is essential to keep all chemicals from contact with the skin, even if they are thought to be harmless. Therefore, use gloves which are rated with a protection index of 6 or with an excellent degradation resistance rating. In all other cases, the gloves should be used for splash protection or short contact only. For more details regarding chemical suitability, please contact Ansell Technical Department. Ensure the chemicals cannot enter via the cuff. 4. Contaminated gloves should be cleaned or washed before removal. 5. The gloves should not come in contact with a naked flame. 6. Gloves shall not be used for protection against ionising radiation nor for use in containment enclosures. 7. GAMMEX® Powder-Free with AMT is intended to be worn after proper surgical disinfection protocol. 8. It should be used by qualified healthcare personnel with sterile procedure training. 9. The glove wearer may be contaminated despite the use of chlorhexidine. 10. GAMMEX® Powder-Free with AMT is not proven to protect against blood borne infection following a needle stick or other skin puncture incident. 11. In case of a glove breach or hole, the user should change the gloves immediately. 12. The gloves should not be used in applications requiring thermal protection. Please ensure the

gloves are used only for the designated purposes. 13. Contents of the pouch are sterile unless the pouch is opened or damaged. Do not use if the pouch is opened or damaged. 14. For Medical Use - Lifetime of the Glove - For normal use, following routine inspection, Ansell recommends Surgical gloves be changed every 60-90 minutes. For gloves used with chemotherapy, we recommend they are changed every 30 minutes, or immediately if damaged.

**INGREDIENTS / HAZARDOUS INGREDIENTS. WARNING:** GAMMEX® Powder-Free with AMT is a sterile, powder-free natural rubber latex surgical glove which incorporates an antimicrobial coating on the inner surface. This antimicrobial coating contains chlorhexidine gluconate (CHG) in addition to zinc salts and other ingredients safe on the skin. For more information, please contact Ansell. GAMMEX® Powder-Free with AMT should not be used by or on anyone sensitized to chlorhexidine. This product also contains natural rubber latex which may cause allergic reactions, including anaphylactic responses. Some gloves might contain other ingredients which are known to be a possible cause of allergies in sensitized persons, who may develop irritant and/or allergic contact reactions. In case the patients experience an anaphylactic or allergic reaction, obtain medical advice immediately.

**PROPER DONNING & DOFFING; How to don gloves:** 1. Remove one glove from the package and inspect it to be sure no pinholes or tears are present. 2. If gloves are ambidextrous, they can be worn on either hand. If not, align the glove's fingers and thumb with the proper hand before donning. 3. Insert five fingers into the cuff and pull the cuff over the wrist. 4. Check for a secure fit around the fingers and palm. The cuff should fit snugly around the wrist. **How to doff gloves:** 1. Grasp the outside edge of the glove near the wrist. 2. Peel the glove away from the hand, turning it inside out. Hold it in the opposite gloved hand. 3. Slide an ungloved finger under the wrist of the remaining glove, being careful not to touch the outside of the glove. 4. Peel the remaining glove off from the inside, creating a "bag" containing both gloves.

**CARE INSTRUCTIONS STORAGE:** Keep away from direct sunlight; store in a dry place at 25°C to a maximum of 40°C and keep in the original packaging. Keep away from ozone sources. If gloves are properly stored, as indicated above, they won't lose their performances and won't change the glove characteristics significantly. If gloves could be affected by ageing or storage, the expiry date is mentioned on the packaging materials. **CLEANING:** This glove is for single-use only and is not designed to be laundered. Re-use of the glove could cause post contamination and postoperative complication. Cleaning and re-sterilization cycle has not been validated for this product by the manufacturer. Product integrity may be compromised by any re-sterilization process used. **DISPOSAL:** After use, the product should be incinerated or disposed of as clinical waste according to the contaminated waste disposal protocol. Used gloves may be contaminated with infectious or other hazardous materials. They should also be disposed once they show any signs of degradation during usage, such as tearing, holes, discoloration and weakening of the gloves.

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