



TRIP-NO-MORE™ CORD COVER

Step-by-Step Guide

For best results, follow these step-by-step instructions when using the Trip-No-More cord cover.



1 Use product within one year of receipt.



2 Floor should be clean and dry before covering cords with the Trip-No-More cord cover.




3 If possible, wait until equipment has been moved into final position for the surgical case before placing Trip-No-More over cords.



4 Apply pressure to outer edges to firmly adhere product to floor.



5 If adhesive strips on underside of the Trip-No-More can't be firmly adhered on both sides of the cords, use a Trip-No-More cord cover with a larger cord channel (see Table 1 below).



6 Make sure the cords being covered are not underneath the adhesive strips on the underside of the product. Cords should be placed within the cord channel as shown and not underneath the adhesive strips.



7 If the associated equipment is being moved during the case, firmly re-adhere cover to the floor by applying pressure to the outer edges, or use a new Trip-No-More cover as needed.



8 Replace the Trip-No-More cover after each surgical case or during the procedure anytime it becomes loose around the edges due to excess moisture.

TRIP-NO-MORE CORD COVER SIZES

The optimal size will depend on the diameter of the cords being covered (select the appropriate channel width) and the total floor area you wish to protect (select the appropriate full cover size). Please see Table 1 for available options. Trip-No-More cord covers can be combined to customize total surface area as needed.

Table 1: Size Comparison Chart

TRIP-NO-MORE CORD FULL COVER SIZE	CORD CHANNEL WIDTH	ITEM NUMBER
8" x 24"	6"	2302
8" x 36"	6"	2303
8" x 125'	6"	2309
8" x 24" (glow in the dark)	6"	2312
8" x 36" (glow in the dark)	6"	2313
10" x 24"	8"	2350
10" x 36"	8"	2351
12" x 24"	10"	2360
12" x 36"	10"	2361
12" x 100'	10"	2369

For product details, instructions for use and an electronic version of this document, visit www.ansell.com/trip-no-more
For education related to all Ansell products, visit www.ansellcares.com

North America: (866) 764-3327 or (818) 534-2500 Email: sandel-cs@ansell.com