

## Sterile, three-layered preparation mat, with leak-proof properties and superb absorption

- **Advanced defenses:** The BioClean™ Chem Prep Mat S-BCPM is composed of three layers – a mesh upper layer allows even liquid distribution, while the middle layer provides excellent absorption, and the impermeable bottom layer protects against leakages
- **Reduced allergy and contamination risks:** This latex-free, low-linting BioClean™ Chem Prep Mat S-BCPM is sterilized by gamma irradiation and individually wrapped to minimize contamination risks
- **Assured quality & optimized cleanliness:** These preparation mats are made in a facility meeting ISO 9001 standards, for product quality you can depend on. They are also packed in an ISO Class 4 cleanroom, ensuring low particulate levels.



### KEY FEATURES & BENEFITS

- **Three-layered construction:** For highly leak-proof absorbency
- **Gamma irradiation sterilization:** Suitable for use in an ISO 4 cleanroom
- **Manufactured in an ISO 9001 facility:** For product quality assurance

### Industries

- Controlled and Critical Environments
- Laboratory and Research
- Production and Manufacturing

## TECHNICAL DATA SHEET

Product Information	
Material	Top layer - Polyolefin apertured film, Middle layer - Air-laid paper, Bottom layer - Polyethylene
Color	Blue
Manufacturing/QMS Audit Standards	Manufacturing QMS Audit Standards ISO 9001
Regulatory	ISO 9001:2015
Packaging	S-BCPM4025B: One piece per sealed inner PE bag; one inner bag per second PE bag; 50 double bags per sealed outer PE bag; four outer bags per lined carton (200 pieces) S-BCPM5640B: One piece per sealed inner PE bag; one inner bag per second PE bag; 25 double bags per sealed outer PE bag; four outer bags per lined carton (100 pieces)
Storage	Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.
Available sizes	40cm x 25cm, 56cm x 40cm
Sterilization Method	GAMMA irradiation (25 kGy)
Sterilization Minimum Dose	25kGy
Sterility Assurance Level	10 <sup>-6</sup>
Cleanroom Class	Class 10/ISO Class 4 & EU GMP Grade A/B and other sterile cleanrooms
Shelf Life	Three (3) years from date of manufacture.
Length	56cm (22") or 40cm (16")
Weight	210gsm

## PHYSICAL / GENERAL PROPERTIES

PROPERTY MEASURED	METHOD USED	UNIT OF MEASUREMENT	TYPICAL VALUE
Sorptive Capacity	IEST-RP-CC004.3	ml/m <sup>2</sup>	852
Sorptive Capacity	IEST-RP-CC004.3	ml/prep mat	83.3



# S-BCPM

Sterile Prep-Mat

## ORDERING INFORMATION

SIZE	40cm x 25cm	56cm x 40cm
REORDER NO.	S-BCPM-4025B	S-BCPM-5640B

## PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE

For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

### Europe, Middle East & Africa Region

Ansell Healthcare Europe NV  
T: +32 (0) 2 528 74 00  
F: +32 (0) 2 528 74 01

### Asia Pacific Region

Ansell Global Trading Center  
T: +603 8310 6688  
F: +603 8310 6699

### North America Region

Ansell Healthcare Products LLC  
US T: +1 800 800 0444  
US F: +1 800 800 0445  
CA T: +1-800-363-8340

### Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.  
T: +52 442 248 1544 / 248 3133

### Australia

Ansell Limited  
T: +61 1800 337 041  
F: +61 1800 803 578

### UK

Ansell Nitritex  
T: +44 1638 663338  
F: +44 1638 668890

## Technology



Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending; [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2022 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

