

## Nitrile elbow-length cleanroom gloves, resistant against a range of chemicals

- **Specialized protection:** BioClean™ Nerva BNAL cleanroom gloves are designed to protect wearers against a range of chemicals, boasting PPE Category III certification
- **Enhanced features:** Their 400mm/16” length enables these BioClean™ gloves to provide elbow-length protection from cross contamination, while a textured surface enhances grip and a beaded cuff ensures sturdier protection with fewer liquid roll-off risks
- **Reduced allergy risk:** The gloves are powder-free, latex-free and accelerator free minimizing allergy risks, notably Type I latex and Type IV chemical allergies and sensitivities



### KEY FEATURES & BENEFITS

- **Chemical resistance protection\*:** Enhancing user protection against accidental chemical splashes
- **Elbow-length design:** Cross contamination protection with extra coverage
- **Non-particulating EasyTearpackaging:** Reduced contamination risks

\*PPE Category III certified

### Industries

- Biotechnology Manufacturing
- Controlled and Critical Environments
- Laboratory Research and Development
- Medical Device Manufacturing
- Pharmaceutical Manufacturing

### Recommended For

- Chemical Splash
- Chemical Handling
- Transferring liquids and solids
- Sample taking and processing



# Nerva BNAL

Non-sterile Disposable Nitrile  
Cleanroom Glove

## TECHNICAL DATA SHEET

Product Information	
Material	Nitrile
Color	White
Shape	Ambidextrous
Cuff	Beaded
Manufacturing/QMS Audit Standards	Manufacturing QMS Audit Standards ISO 9001, NEBB Certified Cleanrooms
Regulatory/Standard Compliance	CE 0493, EN 421:2010, EN 455 Part 2, EN ISO 21420:2020, EN ISO 374-1:2016 Type B (KPT), EN ISO 374-5:2016, Food Contact, Category III, UKCA
Packaging	100 pieces per sealed inner PE bag; one inner PE bag per sealed outer PE bag; 10 outer bags per lined carton (1000 pieces). <b>More sustainable packaging:</b> These gloves are packed in recyclable plastic packaging* and delivered in recycled cardboard shipper cases. *Inner wrap, pouch, bag and liner are made from polyethylene (PE) based film. Always check your local recyclable status as these materials may not be considered suitable for recycling in your location.
Storage	Keep away from direct sunlight; store in a dry place and keep in the original packaging. Keep away from ozone sources. If products are properly stored, as indicated, they won't lose their performances or change characteristics significantly. If products could be affected by ageing or storage, the expiry date is mentioned on the packaging materials.
Country of Origin	Malaysia
Available sizes	S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10), 2XL (10.5-11)
Powder Content	Powder-Free
External Glove Surface	Textured Fingers and Palm
Internal Glove Surface	Chlorinated
Cleanroom Class	Class 10/ISO 4
Shelf Life	Five (5) years from date of manufacture.
Tested for use with Chemotherapy Drugs	No
Protein Level	N/A: contains no natural rubber latex
Anti-static	No



# Nerva BNAL

## Non-sterile Disposable Nitrile Cleanroom Glove

Physical Properties						Testing Method
Sizes	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)	2XL (10.5-11)	EN ISO 21420
Typical Length (mm/in)	400 / 16					
Palm Width (mm/in)	86/3.38	97/3.82	107/4.21	118/4.65	127/5	
Freedom from Holes	1.5 AQL Performance Level 2					EN 374-2
Typical Particle Count $\geq 0.5\mu\text{m}$ (counts / $\text{cm}^2$ )	<1200					IENT-RP-CC005.4
Target Single Wall Palm Thickness (mm/mil)	0.12 / 4.72					EN 455-2
Target Single Wall Finger Thickness (mm/mil)	0.15 / 5.90					EN 455-2
Target Single Wall Cuff Thickness (mm/mil)	0.07 / 2.76					EN 455-2
Force at Break (N) During Aging	$\geq 6$ N					EN 455-2

## IONIC CONTENT

Concentration in $\mu\text{g}/\text{cm}^2$	Typical	Concentration in $\mu\text{g}/\text{cm}^2$	Typical
Ammonium	Available on File	Nitrate	<0.30
Bromide	<0.01	Nitrite	Available on File
Calcium	<0.60	Phosphate	<0.01
Chloride	<0.60	Potassium	<0.10
Fluoride	<0.01	Sodium	<0.10
Lithium	Available on File	Sulphate	<0.10
Magnesium	<0.01	Zinc	<0.10

## ORDERING INFORMATION

SIZE	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)	2XL (10.5-11)
REORDER NO.	BNAL-S	BNAL-M	BNAL-L	BNAL-XL	BNAL-2XL

## PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE



For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

### Europe, Middle East & Africa Region

Ansell Healthcare Europe NV  
T: +32 (0) 2 528 74 00  
F: +32 (0) 2 528 74 01

### Asia Pacific Region

Ansell Global Trading Center  
T: +603 8310 6688  
F: +603 8310 6699

### North America Region

Ansell Healthcare Products LLC  
US T: +1 800 800 0444  
US F: +1 800 800 0445  
CA T: +1-800-363-8340

### Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.  
T: +52 442 248 1544 / 248 3133

### Australia

Ansell Limited  
T: +61 1800 337 041  
F: +61 1800 803 578

### UK

Ansell Nitritex  
T: +44 1638 663338  
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

