

GAMMEX® Latex Moisturizing

S U R G I C A L G L O V E S

Latex powder-free surgical gloves, featuring a built-in moisturizer to rehydrate the skin

- **Soft, supple skin:** GAMMEX® Latex Moisturizing surgical gloves feature a proprietary skin moisturizing technology, enabling a water-based coating with glycerol and dimethicone to retain moisture and rehydrate dry skin, improving the wearer's skin condition
- **Exceptional comfort:** These gloves are also designed using anatomic formers for added comfort and fit
- **Enhanced grip:** With a smooth finish and micro-textured external surface around the palm and fingers, these surgical gloves deliver excellent grip under dry or damp conditions
- **Easy donning:** Ansell's donning aid technology facilitates an easy first don, whether the wearer's hands are dry or damp
- **Optimal allergy prevention:** Additional cleaning processes ensure low residues of latex proteins and chemical accelerators, minimizing the risk of developing latex Type I and chemical Type IV allergies and sensitivities



CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established – this product contains natural rubber latex, which may cause allergic reactions

Recommended For

- Damp donning
- Extra comfort during long, sustained procedures
- Use as an underglove when double gloving

FEATURED TECHNOLOGIES



GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.



GAMMEX[®] Latex Moisturizing

SURGICAL GLOVES

PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder-Free
Color	Natural
Thickness	Standard
Cuff Style	Straight with SUREFIT™ Technology
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Coated with HYDRASOFT™ Technology
Allergy Prevention	None
Grip Level	Moderate
Double Gloving Recommendation	Underglove
Tested For Use With Chemotherapy Drugs	No, not tested in accordance with ASTM D6978
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	AS/NZS 4179, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282, JIS T9107
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Xanthogen Polysulphides, Sodium di-n-butyldithiocarbamate (SDBC), Zinc Diethyldithiocarbamate (ZDEC)

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	330049055	330049060	330049065	330049070	330049075	330049080	330049085	330049090

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 844 494 7854
Email: insidesalesus@ansell.com
Sandel US Tel: 866 764 3327
Sandel email: sandel-cs@ansell.com

Central & South America
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil
Tel: +55-11-3356-3100
Email: latam.medical@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu
Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

India
Tel: +91 22 38124500
Email: info.india@ansell.com
Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com
China
Tel: +86(21) 3827 5005
Email: infochina@ansell.com

ansell.com

Ansell, [®] and [™] are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking.
© 2025 Ansell Limited. All Rights Reserved.