

MICRO-TOUCH® Micro-Thin™ Nitrile

Standard Thin

Nitrile, powder-free medical exam glove with a thin film formulation for excellent tactile sensitivity

- **Heightened sensitivity:** This non-sterile glove is made with a thinner formulation to provide excellent tactile sensitivity and dexterity, suitable for a variety of clinical procedures, including palpating veins and starting IVs
- **Comfortable fit:** This nitrile glove also provides exceptional comfort, fit and feel
- **Secure grip:** The textured surface of the glove promotes sure handling of materials and instruments
- **Smart packaging:** 300 gloves per box reduces restocking and reorders
- **Specialized protection:** This glove has been tested against chemotherapy drugs* and is food handling-compliant

*see product packaging for details

Recommended For

- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Protection from Type I latex allergy in HCW's or patients
- Brief, low-risk examination procedures
- Environmentally friendly facilities

FEATURED TECHNOLOGIES



MICRO-TOUCH® Everyday, Everywhere Protection

MICRO-TOUCH® premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH® Micro-Thin™ Nitrile

Standard Thin

PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Blue	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Finger Thickness (mm/mil)	0.10 / 3.9
Palm Thickness (mm/mil)	0.06 / 2.4	Allergy Prevention	Latex (Type I)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	230 / 9.1	Product Segmentation	Standard Thin

ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	6034310	6034311	6034312	6034313	6034314

PACKAGING AND STORAGE

Packaging	300 gloves per dispenser 10 dispensers per case 3000 gloves per case
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

ASTM D6319, CE 0493, Medical Device Regulation (EU) 2017/745 Class I, FDA21 CFR 177-2600-US Food Contact Approved, Food Type: Rubber, ISO 9001, ISO 13485

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com