

MICRO-TOUCH® DENTA-GLOVE® NRL

D E N T A - G L O V E ® S e r i e s

Latex, powder-free medical exam glove with maximum comfort for a wide range of dental procedures

- **Heightened sensitivity:** The latex formulation and optimized film thickness of this non-sterile glove offers a high level of elasticity, flexibility and sensitivity
- **Easy donning:** The inner chlorinated surface ensures fast and easy donning
- **Suitable for patient contact:** Tested Safe for Oral Examination, ensuring that the glove is safe for patient contact — an ideal choice for dental specialists who switch gloves frequently and experience hand fatigue from continuous use



CAUTION: This product contains natural rubber latex, which may cause allergic reactions

Recommended For

- Enhanced grip in wet environments
- Standard, moderate risk examination procedures
- Oral examinations



MICRO-TOUCH® Everyday, Everywhere Protection

MICRO-TOUCH® premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH® DENTA-GLOVE® NRL

D E N T A - G L O V E ® S e r i e s

PRODUCT INFORMATION

Material	Natural Rubber Latex	Not Made From Natural Rubber Latex	No
Color	Natural, White	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Fully Textured
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.14 / 5.5
Finger Thickness (mm/mil)	0.16 / 6.3	Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)
Tested For Use With Chemotherapy Drugs	No	Sterile	No
Antistatic	Not Tested	Glove Length (mm/inches)	240 / 9.5
Product Segmentation	DENTA-GLOVE® Series		

ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	04653	04654	04655	04656	04657

PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case *90 gloves per dispenser (Size XL) *900 gloves per case (Size XL)
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

ASTM D3578, Category III, EN 420:2003 + A1:2009, EN 455 Part 1, EN 455 Part 2, EN 455 Part 3, EN 455 Part 4, ISO 11193-1:2008, ISO 9001, ISO 13485, ISO 10993-10:2010

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com