

MICRO-TOUCH® SensiClean® II

E s s e n t i a l S e r i e s

Latex medical exam glove offering superior tactile sensitivity and dexterity

- **Enhanced sensitivity:** Thin latex design provides excellent tactile sensitivity and dexterity
- **Comfortable fit:** Elastic natural rubber latex formulation ensures a comfortable feel during wear
- **Reduced skin irritation:** Low residual protein content to reduce risk of contact dermatitis
- **Ideal for double gloving:** The double-chlorinated inner surface is ideal for effortless double gloving while remaining very gentle on fragile skin
- **Secure grip:** The mildly textured surface allows for secure wet and dry grip



CAUTION: This product contains natural rubber latex, which may cause allergic reactions

Recommended For

- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Brief, low-risk examination procedures
- Damp donning



MICRO-TOUCH® Everyday, Everywhere Protection

MICRO-TOUCH® premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH® SensiClean® II

Essential Series

PRODUCT INFORMATION

Material	Natural Rubber Latex	Not Made From Natural Rubber Latex	No
Color	White	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Fully Textured
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.14 / 5.5
Finger Thickness (mm/mil)	0.15 / 5.9	Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)
Tested For Use With Chemotherapy Drugs	No	Sterile	No
Antistatic	Not Tested	Glove Length (mm/inches)	240 / 9.5
Product Segmentation	Essential Series		

ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	4566	4560	4562	4564	4568

PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

EN 455 Part 2, EN 455 Part 3, EN 455 Part 1, EN 420:2003 + A1:2009, ISO 11193, ISO 9001, ISO 13485, ISO 10993-10

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com