

## Clean & non-sterile nitrile RABS/isolator gloves with ultra-low particle level

**Port size:** 305mm/12"

- **Specialized protection:** These nitrile RABS/isolator gloves are autoclave compatible at 121°C for 30 minutes (15 cycles according to internal testing results), resistant to vaporized hydrogen peroxide (VHP) & isopropyl alcohol (IPA)
- **Reduced contamination risk & enhance safety:** These gloves are 100% water leak tested, minimizing risk of contamination
- **Optimized cleanliness:** These cleanroom gloves are processed in accordance with ISO cleanroom standards, for ultra-clean RABS gloves with low particle levels

*For cleaning of isolators with other disinfectants or disinfection method, users should conduct evaluation at their discretion for compatibility or seek Ansell consultation.*



### KEY FEATURES & BENEFITS

- **ASTM D6978-05 compliance:** Suitable for chemotherapy drug handling
- **100% water leak tested:** Assured glove integrity
- **ISO cleanroom processed:** For ultra-clean isolator gloves

### Industries

- Pharmaceutical Manufacturing
- Laboratory Research and Development
- Raw Material and Preparation
- Biotechnology Manufacturing
- Semiconductor manufacturing



## CGL36NIT59

Clean & Non-sterile RABS/Isolator  
Gloves

### TECHNICAL DATA SHEET

|                                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|----------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Reference                      | CGL36NIT59                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Material                               | Nitrile                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Color                                  | White                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Shape                                  | Ambidextrous                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| External Glove Surface                 | Smooth                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Cuff Style                             | Beaded                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Audit Standards                        | ISO 14001, Manufacturing QMS Audit Standards ISO 9001, PPE Regulation 2016 425 Module D, NEBB Certified Cleanrooms                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Packaging Overview                     | <p>One piece per sealed inner PE bag; one inner bag per sealed second inner PE bag; one second inner bag per sealed outer PE bag; 20 outer bags per lined cardboard box (20 pieces).</p> <p><b>More sustainable packaging:</b> These gloves are packed in recyclable plastic packaging* and delivered in recycled cardboard shipper cases.</p> <p>*Inner wrap, pouch, bag and liner are made from polyethylene (PE) based film. Always check your local recyclable status as these materials may not be considered suitable for recycling in your location.</p> |
| Storage Instructions                   | Keep away from direct sunlight; store in a dry place and keep in the original packaging. Keep away from ozone sources. If products are properly stored, as indicated, they won't lose their performances or change characteristics significantly. If products could be affected by ageing or storage, the expiry date is mentioned on the packaging materials.                                                                                                                                                                                                  |
| Country Of Origin                      | Malaysia                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Available Sizes                        | 9.75                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Temperature Range                      | 5 to 40 °C                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Cleanroom Class                        | Class 10/ISO 4                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Shelf Life                             | Three (3) years from date of manufacture.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Tested For Use With Chemotherapy Drugs | Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Protein Level                          | N/A: contains no natural rubber latex                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Standards                              | ASTM D6978, CE 0493, EN 388:2016 + A1:2018, EN 421:2010, EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018, EN ISO 374-5:2016, Food Contact, UKCA                                                                                                                                                                                                                                                                                                                                                                                                                    |



## CGL36NIT59

Clean & Non-sterile RABS/Isolator  
Gloves

### PHYSICAL PROPERTIES

|                                              | Typical Values | Testing Method |
|----------------------------------------------|----------------|----------------|
| Sizes                                        | 9.75           |                |
| Typical Length (mm/in)                       | 840/33         | EN ISO 21420   |
| Palm Width (mm/in)                           | 129/5.1        | EN ISO 21420   |
| Port Size                                    | 305mm/12"      |                |
| Target Single Wall Palm Thickness (mm/mil)   | 0.50/20        | EN 455-2       |
| Target Single Wall Finger Thickness (mm/mil) | 0.60/24        | EN 455-2       |
| Target Single Wall Cuff Thickness (mm/mil)   | 0.50/20        | EN 455-2       |

### PROTECTION AGAINST MECHANICAL RISKS, EN 388:2016+A1:2018

| Risk                     | Performance level |
|--------------------------|-------------------|
| Abrasion Resistance      | Level 4           |
| Blade Cut Resistance     | Level 1           |
| Tear Resistance          | Level 0           |
| Puncture Resistance      | Level 2           |
| ISO 13997 Cut Resistance | Level X           |

### ORDERING INFORMATION

|            |             |            |
|------------|-------------|------------|
|            | PORT SIZE   | 12"/305mm  |
|            | SIZE        | 9.75       |
| CGL36NIT59 | REORDER NO. | CGL36NIT59 |

For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

#### Europe, Middle East & Africa Region

Ansell Healthcare Europe NV  
T: +32 (0) 2 528 74 00  
F: +32 (0) 2 528 74 01

#### Asia Pacific Region

Ansell Global Trading Center  
T: +603 8310 6688  
F: +603 8310 6699

#### North America Region

Ansell Healthcare Products LLC  
US T: +1 800 800 0444  
US F: +1 800 800 0445  
CA T: +1-800-363-8340

#### Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.  
T: +52 442 248 1544 / 248 3133

#### Australia

Ansell Limited  
T: +61 1800 337 041  
F: +61 1800 803 578

#### UK

Ansell Nitritex  
T: +44 1638 663338  
F: +44 1638 668890

### Technology



### Performance Standards and Regulatory Compliance



Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

