

# MICRO-TOUCH® NITRA-TEX™

M u l t i - P u r p o s e

## Nitrile medical exam glove with strong barrier protection

- **Comfortable fit:** Made from a soft synthetic formulation for a comfortable fit and feel
- **Enhanced durability:** Designed for strength and tear resistance
- **Easy donning:** The chlorinated inner coating facilitates easier donning
- **Chemical splash protection:** These gloves resist a range of chemicals used by healthcare professionals
- **Reduced allergy risk:** Nitrile material reduces risk of latex Type I allergies
- **Specialized protection:** Tested and approved for use with chemotherapy drugs\*

\*see product packaging for details

### Recommended For

- Tested for use with chemotherapy drugs
- Standard, moderate risk examination procedures
- Environmentally friendly facilities
- Protection from Type I latex allergy in HCW's or patients

### FEATURED TECHNOLOGIES



## MICRO-TOUCH® Everyday, Everywhere Protection

MICRO-TOUCH® premium medical examination gloves are scientifically designed to complement every need of medical professionals.

Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH® NITRA-TEX™

Multi-Purpose

## PRODUCT INFORMATION

<b>Material</b>	Nitrile	<b>Not Made From Natural Rubber Latex</b>	Yes
<b>Color</b>	Blue	<b>Cuff Length</b>	Standard
<b>Powder Content</b>	Powder-Free	<b>External Glove Surface</b>	Fully Textured
<b>Freedom from Holes (Inspection level I)</b>	1.5 AQL	<b>Finger Thickness (mm/mil)</b>	0.15 / 5.9
<b>Palm Thickness (mm/mil)</b>	0.12 / 4.7	<b>Allergy Prevention</b>	Latex (Type I)
<b>Available Sizes</b>	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	<b>Tested For Use With Chemotherapy Drugs</b>	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
<b>Sterile</b>	No	<b>Antistatic</b>	Not Tested
<b>Glove Length (mm/inches)</b>	250 / 9.8	<b>Product Segmentation</b>	Multi-Purpose

## ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
<b>Product Code</b>	4479	4470	4472	4474	4476

## PACKAGING AND STORAGE

<b>Packaging</b>	100 gloves per dispenser; 10 dispensers per case; 1000 gloves per case
<b>Shelf Life</b>	5 Years
<b>Storage Instructions</b>	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

## STANDARDS AND CERTIFICATIONS

EN 374-2 and -4, EN ISO 21420:2020, ISO 13485, ISO 10993-10, ISO 9001, ISO 11737

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 1-844-494-7854  
 Email: insidesalesus@ansell.com

**Latin America & Caribbean**  
 Tel: +52(442) 296 20 50  
 Email: cslac@ansell.com  
 Brazil Tel: +55-11-3356-3100  
 Brazil Email: luvas.medicas@ansell.com

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: info@ansell.eu

**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com