

## Comfortable gloves offering tactility and protection across numerous applications

- **Elevated comfort:** Ansell MICROFLEX® 92-134 Versatility disposable safety gloves are made of a soft nitrile formulation, which provides a highly comfortable fit and feel while protecting hands
- **Heightened tactility:** These nitrile disposable gloves' lightweight construction provides augmented tactile sensitivity and fine motor control
- **Optimized grip:** These MICROFLEX® gloves are textured on the fingertips to ensure a stronger grip for the safer handling of tools and materials
- **Genuine versatility:** Designed for a broad range of applications, these gloves reinforce hand protection across a wide range of industries



### Industries

- Healthcare
- Automotive
- Food Processing
- Janitorial / Sanitation
- Manufacturing
- Life Sciences
- Government

### Recommended For

- Standard patient examinations
- Picking parts and accessories
- Food processing and handling
- Catering/Meal preparation
- Disposal of waste (trash, biohazard and linens)
- Light assembly tasks
- Lab work, blending, compounding, filling, cleaning
- Preparing pharmaceutical products
- Laboratory analysis

TECHNICAL DATA SHEET

|                                     | Product Information                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|-------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Material                            | Nitrile                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| External Glove Surface              | Textured Fingers                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Audit Standards                     | ISO 13485, ISO 9001                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Product Certification               | Medical Device Regulation (EU) 2017/745 Class I, Personal Protective Equipment Regulation (EU) 2016/425 Category III risks, Korean Standard, ASTM D6319, CE 0493, IFU, EN 16523-1, EN ISO 374-1:2016 Type B JKT, EN ISO 374-5:2016 Virus, EN 455 1-4, EN ISO 21420:2020, FDA21 CFR 177-2600-US Food Contact Approved, US QSR/FDA 510(k) Medical Examination Grade, EU Food Contact Approved, ISO 13485, ISO 18889:2019, ISO 9001, PMDA Notification Number: 13B3X10307AH0001, Trade Name: Ansell Non Latex Medical Examination Glove, UKCA, Anvisa 81496190020 / CA 48.121, INVIMA 2015DM-0013417 |
| Packaging Overview                  | 100 gloves per dispenser, 10 dispensers per case, 1000 gloves per case   90 gloves per dispenser (Size XXL), 900 gloves per case (Size XXL)                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Storage Instructions                | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Country Of Origin                   | China, Malaysia                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Product Segmentation                | Standard Thin                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Antistatic                          | No                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Silicone Free                       | No                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Vulcanization Chemical Accelerators | <ul style="list-style-type: none"><li>Zinc Dibutyldithiocarbamate (ZDBC)</li><li>Zinc Diethyldithiocarbamate (ZDEC)</li></ul> <p><i>Only a very small number of users may be sensitive to this ingredient(s) and hence may develop irritant and/or allergic contact reactions.</i></p>                                                                                                                                                                                                                                                                                                            |

PHYSICAL PROPERTIES

|                                         | Typical Values |             | Testing Method   |
|-----------------------------------------|----------------|-------------|------------------|
| Length (mm/inches)                      | 240 / 9.5      |             | ASTM D6319       |
| Freedom from Holes (Inspection level I) | 1.5 AQL        |             | ASTM D5151       |
| Palm Thickness (mm/mils)                | 0.07 / 2.8     |             | ASTM D6319       |
| Finger Thickness (mm/mils)              | 0.11 / 4.3     |             | ASTM D6319       |
|                                         | BEFORE AGING   | AFTER AGING |                  |
| Ultimate Tensile Strength (MPa)         | ≥ 14           | ≥ 14        | ASTM D412 & D573 |
| Elongation at Break (%)                 | ≥ 500          | ≥ 400       | ASTM D412 & D573 |
| Force at break (N)                      | ≥ 6            | ≥ 6         | EN 455-2         |

ORDERING INFORMATION

| Size         | XS (5.5 - 6) | S (6.5 - 7) | M (7.5 - 8) | L (8.5 - 9) | XL (9.5 - 10) | XXL (10.5 - 11) |
|--------------|--------------|-------------|-------------|-------------|---------------|-----------------|
| Product Code | 92134060     | 92134070    | 92134080    | 92134090    | 92134100      | 92134110        |

Performance Standards and Regulatory Compliance



For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

|                                                                                                                                                                                                                                                                                                                        |                                                                                                                                                                                                                                                                                                                         |                                                                                                                                              |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Europe, Middle East &amp; Africa Region</b><br>Ansell Healthcare Europe NV<br>T: +32 (0) 2 528 74 00<br>F: +32 (0) 2 528 74 01<br><br><b>Asia Pacific Region</b><br>MY +603 8310 6688<br>apac.medical@ansell.com<br>CN: +86 21 38275000<br>infochina@ansell.com<br>JP: +813 5549 8151<br>info.medical.jp@ansell.com | <b>North America Region</b><br>Ansell Healthcare Products LLC<br>T: +1 800 800 0444<br>F: +1 800 800 0445<br><br><b>Latin America &amp; Caribbean Region</b><br>BR +55 (11) 3356-3100<br>latam.medical@ansell.com<br>MX: +52 442 296 2050<br>latam.medical@ansell.com<br>CO: +57 1 288 3247<br>latam.medical@ansell.com | <b>Australia</b><br>Ansell Limited<br>T: +61 1800 337 041<br>F: +61 1800 803 578<br><br><b>Russia</b><br>Ansell PYC<br>Ten. +7 495 258 13 16 |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

