

GAMMEX[®] Latex Underglove

S U R G I C A L G L O V E S

Green latex under gloves for easy breach detection and moisturizing hand protection when double gloving

- **Optimized skin hydration:** These powder-free latex under gloves feature an odorless, water-based skin moisturizer containing glycerol and dimethicone* to help retain moisture and rehydrate skin throughout the duration of use
- **Added safety:** When contrasted with a cream or white-colored outer glove, these green under gloves offer quick and easy breach detection for glove wearer's added protection
- **Effortless donning and doffing:** The internal polymer coating of the GAMMEX[®] Latex Underglove helps to facilitate effortless dry and damp hand donning, whether it be the first don or during intra-operative changes
- **Elevated comfort:** The GAMMEX[®] Latex Underglove is sterile and made using proprietary anatomic-shaped formers and a soft formulation for ergonomic fit and comfort to help minimize hand fatigue during rigorous procedures

*built-in moisturizing agents

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Damp donning
- Easy identification of a glove breach
- Use as an underglove when double gloving

FEATURED TECHNOLOGIES



GAMMEX[®] Innovation that Protects

Driven by innovation, GAMMEX[®] surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

GAMMEX[®] Latex Underglove

SURGICAL GLOVES

PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder-Free
Color	Green
Thickness	Standard
Cuff Style	Straight with SUREFIT™ Technology
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Coated with HYDRASOFT™ Technology
Allergy Prevention	None
Grip Level	Low
Double Gloving Recommendation	Underglove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	AS/NZS 4179, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001, Korean GMP
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Xanthogen Polysulphides, Sodium di-n-butylthiocarbamate (SDBC), Zinc Diethyldithiocarbamate (ZDEC)

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	330050055	330050060	330050065	330050070	330050075	330050080	330050085	330050090

Contact your Ansell representative for ordering or more information.

North America
 US Tel: 800 952 9916
 CA Tel: 844 494 7854
 Email: insidesalesus@ansell.com
 Sandel US Tel: 866 764 3327
 Sandel email: sandel-cs@ansell.com

Central & South America
 Tel: +52(442) 296 20 50
 Email: cslac@ansell.com
Brazil
 Tel: +55-11-3356-3100
 Email: latam.medical@ansell.com

Europe, Middle East & Africa
 Tel: +32 (0) 2 528 74 00
 Email: info@ansell.eu
Asia Pacific
 Tel: +603 8310 6688
 Email: apac.medical@ansell.com

India
 Tel: +91 22 38124500
 Email: info.india@ansell.com
Australia & New Zealand
 Tel: +61 3 9270 7270
 Email: protection@ap.ansell.com
China
 Tel: +86(21) 3827 5005
 Email: infochina@ansell.com