

MICRO-TOUCH® NITRA-TEX™ E.P.

High Risk

Nitrile, powder-free, extended-cuff medical exam glove with strong barrier protection

- **Comfortable fit:** This non-sterile glove is made from a softer synthetic formulation for a more comfortable fit and feel, and offers high tear resistance for optimal protection
- **Easy donning:** The inner chlorinated surface guarantees easy donning
- **Extended protection:** With its extended cuff, this glove is suitable for scenarios where additional wrist protection is required
- **Chemical resistance:** These gloves have excellent resistance to a broad range of commonly used chemicals
- **Reduced allergy risk:** These nitrile gloves also protect from latex Type I allergies
- **Specialized protection:** Tested and approved for use with chemotherapy drugs*




*see product packaging for details

Recommended For

- Extra protection over the wrist and arm
- Sustained, high-risk examination procedures
- Environmentally friendly facilities
- Tested for use with chemotherapy drugs
- Protection from Type I latex allergy in HCW's or patients

FEATURED TECHNOLOGIES



MICRO-TOUCH® Everyday, Everywhere Protection

MICRO-TOUCH® premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH® NITRA-TEX™ E.P.

High Risk

PRODUCT INFORMATION

| | | | |
|---|--|--|--|
| Material | Nitrile | Not Made From Natural Rubber Latex | Yes |
| Color | Blue | Cuff Length | Extended |
| Powder Content | Powder-Free | External Glove Surface | Textured Fingers |
| Freedom from Holes (Inspection level I) | 1.5 AQL | Finger Thickness (mm/mil) | 0.15 / 5.9 |
| Palm Thickness (mm/mil) | 0.11 / 4.3 | Allergy Prevention | Latex (Type I) |
| Available Sizes | XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10) | Tested For Use With Chemotherapy Drugs | Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k) |
| Silicone Free | Not Tested | Sterile | No |
| Antistatic | Not Tested | Double Gloving Recommendation | No |
| Glove Length (mm/inches) | 305 / 12 | Product Segmentation | High Risk |
| Fentanyl Tested | No | | |

PACKAGING AND STORAGE

| | |
|----------------------|--|
| Shelf Life | 3 Years |
| Storage Instructions | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition. |

STANDARDS AND CERTIFICATIONS

ASTM D6319, ASTM D6978, Medical Device Regulation (EU) 2017/745 Class I, EN 455 1-4, ISO 9001, ISO 10993-10:2002, ISO 10993-23, ISO 10993-5, ISO 11737, ISO 13485, Food Contact

Vulcanization Chemical Accelerators

Zinc Dibutyldithiocarbamate (ZDBC)

Only a very small number of users may be sensitive to this ingredient(s) and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email:
insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email:
luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email:
apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email:
protection@ap.ansell.com