

ENCORE® Latex Ortho

S U R G I C A L G L O V E S

Extra-thick latex surgical gloves, designed to withstand the rigors of surgery using power tools

- **Optimal thickness:** On average, ENCORE® Latex Ortho surgical gloves are 50% thicker than standard Ansell latex surgical gloves, providing outstanding hand protection for rigorous ortho procedures, especially from the use of power tools and where bone fragments are present
- **Enhanced grip:** The heavy-textured external finish of these sterile ortho surgical gloves enhances grip performance, facilitating the secure handling of power tools in both dry and wet conditions
- **Effortless donning:** ENCORE® Latex Ortho gloves are powder-free, with a unique donning aid technology, featuring a polymer-lined inner coating, for speedy and easy dry and damp hand donning, whether it is the first don or during intra-operative changes
- **Chemo-tested for safe use**:** Tested for safe use with chemotherapy drugs in accordance with ASTM D6978 requirements

**contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Enhanced grip in wet environments
- Damp donning

FEATURED TECHNOLOGIES

 **A.R.T.**
Pinhole Barrier Technology

 **DERMASHIELD™**
Donning Aid Technology



ENCORE® The Confidence to Perform

Dedicated to delivering optimized solutions engineered to provide a high degree of confidence to surgeons, O.R. nurses and surgical technicians seeking protection with a high level of dexterity, comfort and grip so they can perform at their best.



ENCORE[®] Latex Ortho

SURGICAL GLOVES

PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder-Free
Color	Brown
Thickness	Ortho
Cuff Style	Beaded with SUREFIT™ Technology
External Glove Surface	Textured Fingers and Palm
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology
Allergy Prevention	None
Grip Level	High
Double Gloving Recommendation	Outer glove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	ASTM D3577, ASTM D6978, EN 16523-1, EN 388:2016 +A1:2018, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN 374-2 and -4, EN ISO 374-5:2016, ISO 10282
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC)

ORDERING INFORMATION

Size	6	6.5	7	7.5	8	8.5	9
Product Code	330106060	330106065	330106070	330106075	330106080	330106085	330106090

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 844 494 7854
Email: insidesalesus@ansell.com
Sandel US Tel: 866 764 3327
Sandel email: sandel-cs@ansell.com

Central & South America
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil
Tel: +55-11-3356-3100
Email: latam.medical@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu
Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

India
Tel: +91 22 38124500
Email: info.india@ansell.com
Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com
China
Tel: +86(21) 3827 5005
Email: infochina@ansell.com

ansell.com

Ansell, [®] and [™] are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking.
© 2025 Ansell Limited. All Rights Reserved.