

Sterile Singles

Synthetic, sterile single packed, powder-free examination glove ideal for brief, sterile applications

Single-packed, synthetic medical exam glove ideal for brief, sterile applications

- **Sterile:** Packaged as a single sterile glove per pouch, this powder-free glove is ideal for sterile medical exam (non-surgical) procedures
- **Reduced allergy risk:** This synthetic glove is made from copolymer materials (ethylene vinyl acetate film), making it suitable for use by individuals with allergies to natural rubber latex (Type I allergies)
- **Double gloving-friendly:** ETHIPARAT Sterile Singles can also be used as glove liners for increased protection
- **Cost-effective:** These synthetic gloves for sterile use are offered at a competitive price point
- **Easy donning:** These gloves are packaged on backing paper, making them easy to don



Industries

- Healthcare

Recommended For

- Sterile procedures outside of the O.R.
- Brief, low-risk examination procedures
- Protection from Type I latex allergy in HCW's or patients

TECHNICAL DATA SHEET

	Product Information
External Glove Surface	Smooth
Audit Standards	ISO 13485, ISO 9001
Product Certification	MDD 93/42/EWG Klasse I Steril
Packaging Overview	100 gloves per dispenser 10 dispensers per case 1000 gloves per case
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Country Of Origin	Taiwan
Product Segmentation	Essential Series
Antistatic	Not Tested
Vulcanization Chemical Accelerators	<ul style="list-style-type: none">None

PHYSICAL PROPERTIES

	Typical Values		Testing Method
Length (mm/inches)	278 / 10.9		ASTM D3767,EN 420
Freedom from Holes (Inspection level I)	1.5 AQL		ASTM D5151,EN 455-1
Palm Thickness (mm/mils)	0.03 / 1.2		ASTM D3767,EN 420
Finger Thickness (mm/mils)	0.03 / 1.2		ASTM D3767,EN 420
	BEFORE AGING	AFTER AGING	
Force at break (N)	≥2	≥2	EN 455-2

Performance Standards and Regulatory Compliance



For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region
Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region
MY +603 8310 6688
apac.medical@ansell.com
CN: +86 21 38275000
infochina@ansell.com
JP: +813 5549 8151
info.medical.jp@ansell.com

North America Region
Ansell Healthcare Products LLC
T: +1 800 800 0444
F: +1 800 800 0445

Latin America & Caribbean Region
BR +55 (11) 3356-3100
latam.medical@ansell.com
MX: +52 442 296 2050
latam.medical@ansell.com
CO: +57 1 288 3247
latam.medical@ansell.com

Australia
Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

Russia
Ansell PYC
Ten. +7 495 258 13 16

Technology



Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking © 2024 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

