

GAMMEX[®] Latex Underglove

S U R G I C A L G L O V E S

Green latex under gloves for easy breach detection and moisturizing hand protection when double gloving

- **Optimized skin hydration:** These powder-free latex under gloves feature an odorless, water-based skin moisturizer containing glycerol and dimethicone* to help retain moisture and rehydrate skin throughout the duration of use
- **Added safety:** When contrasted with a cream or white-colored outer glove, these green under gloves offer quick and easy breach detection for glove wearer's added protection
- **Effortless donning and doffing:** The internal polymer coating of the GAMMEX[®] Latex Underglove helps to facilitate effortless dry and damp hand donning, whether it be the first don or during intra-operative changes
- **Elevated comfort:** The GAMMEX[®] Latex Underglove is sterile and made using proprietary anatomic-shaped formers and a soft formulation for ergonomic fit and comfort to help minimize hand fatigue during rigorous procedures

*built-in moisturizing agents

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Damp donning
- Easy identification of a glove breach
- Use as an underglove when double gloving

FEATURED TECHNOLOGIES



GAMMEX[®] Innovation that Protects

Driven by innovation, GAMMEX[®] surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

GAMMEX[®] Latex Underglove

SURGICAL GLOVES

PRODUCT INFORMATION

| | |
|--|--|
| Material | Natural Rubber Latex |
| Powder Content | Powder-Free |
| Color | Green |
| Thickness | Standard |
| Cuff Style | Straight with SUREFIT™ Technology |
| External Glove Surface | Smooth with Micro-Textured Finish |
| Internal Glove Surface | Coated with HYDRASOFT™ Technology |
| Allergy Prevention | None |
| Grip Level | Low |
| Double Gloving Recommendation | Underglove |
| Tested For Use With Chemotherapy Drugs | Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k) |
| Freedom from Holes (Inspection level I) | 0.65 AQL |
| Sterile | Yes |
| Product Standards | AS/NZS 4179, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282 |
| Quality / Environmental Standards | EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001, Korean GMP |
| Packaging | 50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case |
| Shelf Life | 3 years |
| Storage Instructions | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition. |
| Vulcanization Chemical Accelerators | Xanthogen Polysulphides, Sodium di-n-butylthiocarbamate (SDBC), Zinc Diethyldithiocarbamate (ZDEC) |

ORDERING INFORMATION

| Size | 5.5 | 6 | 6.5 | 7 | 7.5 | 8 | 8.5 | 9 |
|---------------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| Product Code | 330050055 | 330050060 | 330050065 | 330050070 | 330050075 | 330050080 | 330050085 | 330050090 |

Contact your Ansell representative for ordering or more information.

North America
 US Tel: 800 952 9916
 CA Tel: 844 494 7854
 Email: insidesalesus@ansell.com
 Sandel US Tel: 866 764 3327
 Sandel email: sandel-cs@ansell.com

Central & South America
 Tel: +52(442) 296 20 50
 Email: cslac@ansell.com
Brazil
 Tel: +55-11-3356-3100
 Email: latam.medical@ansell.com

Europe, Middle East & Africa
 Tel: +32 (0) 2 528 74 00
 Email: info@ansell.eu
Asia Pacific
 Tel: +603 8310 6688
 Email: apac.medical@ansell.com

India
 Tel: +91 22 38124500
 Email: info.india@ansell.com
Australia & New Zealand
 Tel: +61 3 9270 7270
 Email: protection@ap.ansell.com
China
 Tel: +86(21) 3827 5005
 Email: infochina@ansell.com