

# MEDI-GRIP® Latex

S U R G I C A L   G L O V E S

## Latex, powder-free surgical glove with a unique external surface treatment for a safe and secure grip

### KEY FEATURES & BENEFITS

- Good handling and grip achieved through a unique external surface treatment
- Soft and malleable latex formulation provides ample dexterity and good comfort
- Made using anatomic formers to enhance comfort and fit

*CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.*

### Recommended For

- Brief, low-risk examination procedures
- Extra comfort during long, sustained procedures
- Sterile procedures in the O.R.
- Use as an outer glove when double gloving



**MEDI-GRIP®**  
Everyday Protection you can  
Depend on

Trusted around the world to deliver the every day protection surgeons, O.R. nurses and surgical technicians can depend on to perform their daily tasks safely and comfortably.

# MEDI-GRIP® Latex

## SURGICAL GLOVES

### PRODUCT INFORMATION

|                                         |                                                                                                                             |
|-----------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Material                                | Natural Rubber Latex                                                                                                        |
| Powder Content                          | Powder-Free                                                                                                                 |
| Color                                   | Natural                                                                                                                     |
| Thickness                               | Standard                                                                                                                    |
| Cuff Style                              | Beaded                                                                                                                      |
| External Glove Surface                  | Smooth                                                                                                                      |
| Internal Glove Surface                  | Polymer Coated                                                                                                              |
| Allergy Prevention                      | None                                                                                                                        |
| Grip Level                              | Moderate                                                                                                                    |
| Double Gloving Recommendation           | Outer glove                                                                                                                 |
| Tested For Use With Chemotherapy Drugs  | No, not tested in accordance with ASTM D6978                                                                                |
| Freedom from Holes (Inspection level I) | 0.65 AQL                                                                                                                    |
| Sterile                                 | Yes                                                                                                                         |
| Product Standards                       | ASTM D3577, EN 16523-1, EN 374-2 and -4, EN 420:2003 + A1:2009, EN 455 1-4, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282 |
| Quality / Environmental Standards       | EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001                                                                    |
| Packaging                               | 50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case                                                        |
| Shelf Life                              | 3 years                                                                                                                     |
| Storage Instructions                    | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.                    |
| Vulcanization Chemical Accelerators     | Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC)                                                      |

### ORDERING INFORMATION

| Size         | 5.5       | 6         | 6.5       | 7         | 7.5       | 8         | 8.5       | 9         |
|--------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| Product Code | 330200955 | 330200960 | 330200965 | 330200970 | 330200975 | 330200980 | 330200985 | 330200990 |

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: [insidesalesus@ansell.com](mailto:insidesalesus@ansell.com)  
 Sandel US Tel: 866 764 3327  
 Sandel email: [sandel-cs@ansell.com](mailto:sandel-cs@ansell.com)

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: [cslac@ansell.com](mailto:cslac@ansell.com)  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: [latam.medical@ansell.com](mailto:latam.medical@ansell.com)

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: [info@ansell.eu](mailto:info@ansell.eu)  
**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: [apac.medical@ansell.com](mailto:apac.medical@ansell.com)

**India**  
 Tel: +91 22 38124500  
 Email: [info.india@ansell.com](mailto:info.india@ansell.com)  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: [protection@ap.ansell.com](mailto:protection@ap.ansell.com)  
**China**  
 Tel: +86(21) 3827 5005  
 Email: [infochina@ansell.com](mailto:infochina@ansell.com)