

# MICRO-TOUCH® Blue Nitrile

Standard Thin

## Nitrile, powder-free medical exam glove with exceptional fit and reliable protection

- **Comfortable fit:** These non-sterile gloves features an advanced stretch formulation for ultimate comfort, fit and feel
- **Heightened sensitivity:** Advanced formulation provides better tactile sensitivity
- **Strong barrier integrity:** An AQL better than required by industry standard ASTM D6319 delivers consistent, high-quality production
- **Smartly packed:** 200-count packaging reduces waste
- **Specialized protection:** Tested for use with 17 chemicals and approved for use with chemotherapy drugs, making it ideal for nearly all clinical settings




\*see product packaging for details

### Recommended For

- Tested for use with chemotherapy drugs
- Protection from Type I latex allergy in HCW's or patients
- Standard, moderate risk examination procedures
- Environmentally friendly facilities

### FEATURED TECHNOLOGIES



## MICRO-TOUCH® Everyday, Everywhere Protection

MICRO-TOUCH® premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH® Blue Nitrile

Standard Thin

## PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Blue	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Finger Thickness (mm/mil)	0.09 / 3.5
Palm Thickness (mm/mil)	0.07 / 2.8	Allergy Prevention	Latex (Type I)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	240 / 9.5	Product Segmentation	Standard Thin

## ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	313041060	313041065	313041070	313041075	313041080

## PACKAGING AND STORAGE

Packaging	200 gloves per dispenser 10 dispensers per case 2000 gloves per case
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

## STANDARDS AND CERTIFICATIONS

ASTM D6319, CE 0493, Category III, Medical Device Regulation (EU) 2017/745 Class I, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, FDA21 CFR 177-2600-US Food Contact Approved, ISO 10993-10:2002, ISO 9001, ISO 13485, Food Contact, UKCA

## Vulcanization Chemical Accelerators

Zinc Dibutyldithiocarbamate (ZDBC)	Zinc Diethyldithiocarbamate (ZDEC)
------------------------------------	------------------------------------

Only a very small number of users may be sensitive to this ingredient(s) and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email: insidesalesus@ansell.com

Latin America & Caribbean  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

Asia Pacific  
Tel: +603 8310 6688  
Email: apac.medical@ansell.com

Australia & New Zealand  
Tel: +61 3 9270 7270  
Email: protection@ap.ansell.com