

MICRO-TOUCH® NITRAFREE™

High Touch

Provides comfortable and reliable hand protection

- **Reduced allergy risk:** Manufactured without accelerators, helping to protect users from Type IV chemical allergies
- **Comfortable fit:** These nitrile medical exam gloves are made from a durable and soft nitrile formulation for enhanced comfort
- **Heightened tactility:** Thin and lightweight to ensure excellent tactility while working
- **Specialized protection:** Tested and approved for use with chemotherapy drugs*

*see product packaging for details



Recommended For

- Tested for use with chemotherapy drugs
- Standard patient examinations
- Protection from latex Type I allergy and minimize chemical Type IV allergies in HCW's or patients

FEATURED TECHNOLOGIES



MICRO-TOUCH® Everyday, Everywhere Protection

MICRO-TOUCH® premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH® NITRAFREE™

High Touch

PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Pink	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Finger Thickness (mm/mil)	0.12 / 5.5
Palm Thickness (mm/mil)	0.08 / 3.1	Allergy Prevention	Latex (Type I), Chemical (Type IV)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Silicone Free	Not Tested	Sterile	No
Antistatic	Not Tested	Double Gloving Recommendation	Yes, as underglove
Glove Length (mm/inches)	240 / 9.5	Product Segmentation	High Touch

ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	6034510	6034511	6034512	6034513	6034514

PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case
Shelf Life	5 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

ASTM D6319, Category III, CE 0493, Medical Device Regulation (EU) 2017/745 Class I, EN 16523-1, EN 421:2010, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, FDA21 CFR 177-2600-US Food Contact Approved, ISO 9001, ISO 13485, ISO 10993-10:2002, UKCA

Vulcanization Chemical Accelerators

Made without accelerators

Only a very small number of users may be sensitive to this ingredient(s) and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com