

MICRO-TOUCH[®] DermaClean[®]

Essential Series

Latex, powder-free medical exam glove with excellent grip properties


- **Secure grip:** Thanks to its textured outer surface, this non-sterile glove has excellent grip levels in wet and dry conditions, and is ideal for healthcare providers who handle instrumentation, test tubes or infusion equipment
- **Comfortable fit:** Ansell's natural rubber latex formulation ensures flexibility in a comfortable, well-fitting glove
- **Easy donning:** The smooth, chlorinated inner surface facilitates ease of donning



CAUTION: This product contains natural rubber latex, which may cause allergic reactions

Recommended For

- Brief, low-risk examination procedures



MICRO-TOUCH[®] Everyday, Everywhere Protection

MICRO-TOUCH[®] premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH® DermaClean®

Essential Series

PRODUCT INFORMATION

Material	Natural Rubber Latex	Not Made From Natural Rubber Latex	No
Color	White	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Fully Textured
Freedom from Holes (Inspection level I)	1.5 AQL	Finger Thickness (mm/mil)	0.16 / 6.3
Palm Thickness (mm/mil)	0.14 / 5.5	Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)
Tested For Use With Chemotherapy Drugs	No	Silicone Free	Yes
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	240 / 9.5	Product Segmentation	Essential Series

ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	4570	4572	4574	4576	4578

PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case *80 gloves per dispenser (Size XL) *800 gloves per case (Size XL)
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

CA Prop 65, CE Number, Medical Device Regulation (EU) 2017/745 Class I, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 13485, ISO 10993-5, ISO 9001, ISO 11193-1:2008, ISO 10993-1, ISO 16604, ISO 10993-10, ISO 10993-10:2002, Food Type: Rubber

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com