

# ENCORE<sup>®</sup> Latex Underglove

S U R G I C A L   G L O V E S

## Powder-free latex undergloves, designed for double gloving during surgical procedures

- **Enhanced safety:** The ENCORE<sup>®</sup> Latex Underglove is surgical green in color, contrasting with most outer gloves and making glove breaches easier to detect
- **Improved comfort and hydration:** These powder-free latex gloves also feature an internal coating, designed with HYDRASOFT<sup>™</sup> Skin Moisturizing Technology\*\*, retaining moisture and rehydrating the skin
- **Effortless donning and doffing:** Their DERMASHIELD<sup>™</sup> Technology inner coating and low-friction outer surface make these sterile gloves easy and comfortable to don and remove in dry and damp conditions, facilitating double gloving and intra-operative changes

\*\*contains glycerol and dimethicone moisturizing agents

**CAUTION:** Safe use of these gloves by latex-sensitized individuals has not been established — this product contains natural rubber latex, which may cause allergic reactions

### Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Damp donning
- Easy identification of a glove breach
- Use as an underglove when double gloving

### FEATURED TECHNOLOGIES



## ENCORE<sup>®</sup> The Confidence to Perform

Dedicated to delivering optimized solutions engineered to provide a high degree of confidence to surgeons, O.R. nurses and surgical technicians seeking protection with a high level of dexterity, comfort and grip so they can perform at their best.



# ENCORE® Latex Underglove

## SURGICAL GLOVES

### PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder-Free
Color	Green
Thickness	Micro
Cuff Style	Beaded with SUREFIT™ Technology
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Coated with HYDRASOFT™ Technology
Allergy Prevention	None
Grip Level	Low
Double Gloving Recommendation	Underglove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	ASTM D3577, ASTM D7160, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN 374-2 and -4, EN ISO 374-5:2016, ISO 10282
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Zinc Dibutyldithiocarbamate (ZDBC), Zinc Diethyldithiocarbamate (ZDEC)

### ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	2018455	2018460	2018465	2018470	2018475	2018480	2018485	2018490

Contact your Ansell representative for ordering or more information.

**North America**  
US Tel: 800 952 9916  
CA Tel: 844 494 7854  
Email: [insidesalesus@ansell.com](mailto:insidesalesus@ansell.com)  
Sandel US Tel: 866 764 3327  
Sandel email: [sandel-cs@ansell.com](mailto:sandel-cs@ansell.com)

**Central & South America**  
Tel: +52(442) 296 20 50  
Email: [cslac@ansell.com](mailto:cslac@ansell.com)  
**Brazil**  
Tel: +55-11-3356-3100  
Email: [latam.medical@ansell.com](mailto:latam.medical@ansell.com)

**Europe, Middle East & Africa**  
Tel: +32 (0) 2 528 74 00  
Email: [info@ansell.eu](mailto:info@ansell.eu)  
**Asia Pacific**  
Tel: +603 8310 6688  
Email: [apac.medical@ansell.com](mailto:apac.medical@ansell.com)

**India**  
Tel: +91 22 38124500  
Email: [info.india@ansell.com](mailto:info.india@ansell.com)  
**Australia & New Zealand**  
Tel: +61 3 9270 7270  
Email: [protection@ap.ansell.com](mailto:protection@ap.ansell.com)  
**China**  
Tel: +86(21) 3827 5005  
Email: [infochina@ansell.com](mailto:infochina@ansell.com)

**ansell.com**

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking).  
© 2025 Ansell Limited. All Rights Reserved.