

ENCORE[®] Perry Style 42[™] PF

S U R G I C A L G L O V E S

Latex powder-free surgical gloves, ideally combining protection and tactile sensitivity

- **Balanced performance:** Thanks to their optimal thickness levels, ENCORE[®] Perry Style 42[™] PF latex surgical gloves ideally combine the glove wearer's need for barrier protection and tactile sensitivity, offering better instrument 'feel' during surgery
- **Improved comfort:** The proprietary anatomic formers used in the design of these ENCORE[®] gloves ensure a comfortable fit and minimal hand fatigue with sustained use
- **Effortless donning and doffing:** These latex surgical gloves' internal surface has a special polymer coating, making them swift and easy to don
- **Chemo-tested:** They are also tested for use with chemotherapy drugs in accordance with ASTM D6978*

*Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

Recommended For

- Brief, low-risk examination procedures
- Sterile procedures in the O.R.
- Use as an outer glove when double gloving
- Sterile procedures outside of the O.R.



ENCORE[®] The Confidence to Perform

Dedicated to delivering optimized solutions engineered to provide a high degree of confidence to surgeons, O.R. nurses and surgical technicians seeking protection with a high level of dexterity, comfort and grip so they can perform at their best.



ENCORE[®] Perry Style 42[™] PF

SURGICAL GLOVES

PRODUCT INFORMATION

| | |
|---|--|
| Material | Natural Rubber Latex |
| Powder Content | Powder-Free |
| Color | Natural |
| Thickness | Standard |
| Cuff Style | Beaded |
| External Glove Surface | Smooth with Micro-Textured Finish |
| Internal Glove Surface | Polymer Coated |
| Allergy Prevention | None |
| Grip Level | Moderate |
| Double Gloving Recommendation | Outer glove |
| Tested For Use With Chemotherapy Drugs | Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k) |
| Freedom from Holes (Inspection level I) | 0.65 AQL |
| Sterile | Yes |
| Product Standards | ASTM D3577, COFERPIS: 1243C2023 SSA, EN 420:2003 + A1:2009 |
| Quality / Environmental Standards | ISO 13485, ISO 9001 |
| Packaging | 50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case |
| Shelf Life | 3 years |
| Storage Instructions | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition. |
| Vulcanization Chemical Accelerators | Zinc Dibutyldithiocarbamate (ZDBC), Zinc Diethyldithiocarbamate (ZDEC) |

ORDERING INFORMATION

| Size | 5.5 | 6 | 6.5 | 7 | 7.5 | 8 | 8.5 | 9 |
|--------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| Product Code | 5711100PF | 5711101PF | 5711102PF | 5711103PF | 5711104PF | 5711105PF | 5711106PF | 5711107PF |

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 844 494 7854
Email: insidesalesus@ansell.com
Sandel US Tel: 866 764 3327
Sandel email: sandel-cs@ansell.com

Central & South America
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil
Tel: +55-11-3356-3100
Email: latam.medical@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu
Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

India
Tel: +91 22 38124500
Email: info.india@ansell.com
Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com
China
Tel: +86(21) 3827 5005
Email: infochina@ansell.com

ansell.com

Ansell, [®] and [™] are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking.
© 2025 Ansell Limited. All Rights Reserved.