

# GAMMEX® Latex Ortho

S U R G I C A L   G L O V E S

## Ultra-thick powder-free latex gloves designed to withstand rigors of surgery using power tools

- **Added protection:** These brown color surgical gloves are powder-free and made 50% thicker\* than Ansell's standard latex gloves using Ansell's extra-thick formulation, making them suitable for high-risk procedures such as ortho surgery to protect wearers from power tools and bone fragments
- **Exceptional grip:** GAMMEX® Latex Ortho powder-free sterile gloves feature textured fingers and palm for secure grip when handling power tools under dry or damp conditions
- **Effortless donning:** Ansell's unique donning aid technology enables a polymer coating to be applied on the inside of the glove for easy donning and doffing, be it the first don or during intraoperative changes
- **Chemo-tested:** These sterile latex surgical gloves are tested per ASTM D6978 and may be used with a wide range of cytotoxic agents\*\*

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

*\*Compared to Ansell's standard latex gloves*

*\*\*Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations*

### Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Damp donning
- Enhanced grip in wet environments

### FEATURED TECHNOLOGIES



## GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX<sup>®</sup> Latex Ortho

SURGICAL GLOVES

## PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder-Free
Color	Brown
Thickness	Ortho
Cuff Style	Beaded with SUREFIT™ Technology
External Glove Surface	Textured Fingers and Palm
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology
Allergy Prevention	None
Grip Level	High
Double Gloving Recommendation	Outer glove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	AS/NZS 4179, ASTM D6978, EN 16523-1, EN 388:2016 +A1:2018, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001, Korean GMP
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC)

## ORDERING INFORMATION

Size	6	6.5	7	7.5	8	8.5	9
Product Code	330065060	330065065	330065070	330065075	330065080	330065085	330065090

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: insidesalesus@ansell.com  
 Sandel US Tel: 866 764 3327  
 Sandel email: sandel-cs@ansell.com

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: cslac@ansell.com  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: [latam.medical@ansell.com](mailto:latam.medical@ansell.com)

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: info@ansell.eu  
**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**India**  
 Tel: +91 22 38124500  
 Email: info.india@ansell.com  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com  
**China**  
 Tel: +86(21) 3827 5005  
 Email: infochina@ansell.com