

Isolator glove assembly unit, designed for fixing gloves to isolator sleeves

- **Enhanced convenience and reliability:** Made of PVC and silicone, the BioClean™ Isolator Sleeve/Glove Assembly is specially designed to enable the easy connection of gloves to isolator sleeves
- **Improved hygiene:** This isolator glove assembly unit is double-bagged in EasyTear packaging, ensuring hygiene and convenience

KEY FEATURES & BENEFITS

- **PVC C ring and silicone ring:** For easy connection of gloves to isolator sleeves
- **Double-bagged in EasyTear packaging:** Clean and convenient





BBCO 100

Isolator Glove Assembly Unit

TECHNICAL DATA SHEET

Product Reference	BBCO 100
Material	PVC, Silicone
Color	Grey, Red
Packaging Overview	One assembly unit consisting of 1 x C ring and 1 x O ring per inner PE bag; one inner bag per outer bag, five outer bags per shipper case.
Storage Instructions	Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.
Country Of Origin	Vietnam
Available Sizes	11
Cleanroom Class	Class 10/ISO 4
Shelf Life	Five (5) years from date of manufacture.



BBCO 100

Isolator Glove Assembly Unit

ORDERING INFORMATION

	PORT SIZE	
	SIZE	11
BBCO 100	REORDER NO.	BBCO-100

For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
US T: +1 800 800 0444
US F: +1 800 800 0445
CA T: +1-800-363-8340

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

UK

Ansell Nitritex
T: +44 1638 663338
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending; www.ansell.com/patentmarking © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

