

### GAMMEX® Latex Ortho

SURGICAL GLOVES

# Ultra-thick powder-free latex gloves designed to withstand rigors of surgery using power tools

- Added protection: These brown color surgical gloves are powder-free and made 50% thicker\* than Ansell's standard latex gloves using Ansell's extrathick formulation, making them suitable for high-risk procedures such as ortho surgery to protect wearers from power tools and bone fragments
- Exceptional grip: GAMMEX<sup>®</sup> Latex Ortho powder-free sterile gloves feature textured fingers and palm for secure grip when handling power tools under dry or damp conditions
- Effortless donning: Ansell's unique donning aid technology enables a
  polymer coating to be applied on the inside of the glove for easy donning
  and doffing, be it the first don or during intraoperative changes
- Chemo-tested: These sterile latex surgical gloves are tested per ASTM D6978 and may be used with a wide range of cytotoxic agents\*\*

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

\*Compared to Ansell's standard latex gloves

\*\*Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

#### Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Damp donning
- Enhanced grip in wet environments

#### FEATURED TECHNOLOGIES







## GAMMEX<sup>®</sup> Innovation that Protects

Driven by innovation, GAMMEX<sup>®</sup> surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.



## GAMMEX® Latex Ortho

SURGICAL GLOVES

Material	Natural Rubber Latex						
Powder Content	Powder-Free						
Color	Brown						
Thickness	Ortho						
Cuff Style	Beaded with SUREFIT™ Technology						
External Glove Surface	Textured Fingers and Palm						
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology						
Allergy Prevention	None						
Grip Level	High						
Double Gloving Recommendation	Outer glove						
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)						
Freedom from Holes (Inspection level I)	0.65 AQL						
Sterile	Yes						
Product Standards	AS/NZS 4179, ASTM D6978, EN 16523-1, EN 388:2016 +A1:2018, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016						
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001, Korean GMP						
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case						
Shelf Life	3 years						
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.						
Vulcanization Chemical Accelerators	Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC)						
ORDERING INFORMATIO	N						
Size	6	6.5	7	7.5	8	8.5	9
Product Code	330065060	330065065	330065070	330065075	330065080	330065085	33006509

#### Contact your Ansell representative for ordering or more information.

North America US Tel: 800 952 9916 CA Tel: 844 494 7854 Email: insidesalesus@ansell.com Sandel US Tel: 866 764 3327 Sandel email: sandel-cs@ansell.com

Central & South America
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil
Tel: +55-11-3356-3100
Email: latam.medical@ansell.com

Europe, Middle East & Africa Tel: +32 (0) 2 528 74 00 Email: info@ansell.eu Asia Pacific Tel: +603 8310 6688 Email: apac.medical@ansell.com

India
Tet: +91 22 38124500
Email: info.india@ansell.com
Australia & New Zealand
Tet: +61 3 9270 7270
Email: protection@ap.ansell.com
China
Tet: +86(21) 3827 5005
Email: infochina@ansell.com