

GAMMEX® Non-Latex PI Micro

S U R G I C A L G L O V E S

Polyisoprene, skin-friendly surgical gloves, for exceptional dexterity and tactile sensitivity

- **Elevated tactile sensitivity and dexterity:** The combination of Ansell's unique soft, yet strong, polyisoprene formulation and process helps make GAMMEX® Non-Latex PI Micro surgical gloves 15% thinner*, delivering exceptional surgical finesse and hand protection
- **Skin-friendly:** These non-latex surgical gloves eliminate any risk of developing a Type I latex allergy, while their ground-breaking non-sensitizing polyisoprene (PI) formulation eliminates the use of chemical accelerators known to cause allergic contact dermatitis
- **Excellent grip efficiency:** The unique micro-textured external finish offers safe grip when handling surgical instruments under both dry and damp conditions
- **Chemo-tested:** Having been tested per ASTM D6978, these sterile surgical gloves may be safely used with a range of chemotherapy drugs***



***Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

Recommended For

- All surgeries, especially delicate procedures in ophthalmology, micro, and cardiovascular
- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Protection from chemical accelerators such as thiurams, carbamates and mercaptos known to cause chemical Type IV allergies
- Having passed the Mouse Embryo Assay (MEA) test, may be used for In Vitro Fertilization (IVF) procedures

FEATURED TECHNOLOGIES



GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.



GAMMEX® Non-Latex PI Micro

SURGICAL GLOVES

PRODUCT INFORMATION

Material	PI (Polyisoprene)
Powder Content	Powder-Free
Color	White
Thickness	Micro
Cuff Style	Beaded with SUREFIT™ Technology
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology
Allergy Prevention	Latex (Type I)
Grip Level	Moderate
Double Gloving Recommendation	Outer glove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	AS/NZS 4179, ASTM D3577, ASTM D7160, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, JIS T9107
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001, Korean GMP
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Xanthates, Zinc Dialkyldithiocarbamate

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	340057055	340057060	340057065	340057070	340057075	340057080	340057085	340057090

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 844 494 7854
Email: insidesalesus@ansell.com
Sandel US Tel: 866 764 3327
Sandel email: sandel-cs@ansell.com

Central & South America
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil
Tel: +55-11-3356-3100
Email: latam.medical@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu
Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

India
Tel: +91 22 38124500
Email: info.india@ansell.com
Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com
China
Tel: +86(21) 3827 5005
Email: infochina@ansell.com

ansell.com

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking.
© 2025 Ansell Limited. All Rights Reserved.