

### **GAMMEX**® Non-Latex

SURGICAL GLOVES

# Neoprene, chemical accelerator-free surgical gloves offer advanced allergy protection

- Advanced allergy protection: GAMMEX<sup>®</sup> Non-Latex neoprene surgical gloves are synthetic and eliminate risks of Type I latex allergy. These sterile gloves are made 100% free of chemical accelerators and minimize any risk of sensitization to allergens leading to allergic contact dermatitis
- **Proven durability:** Strong neoprene glove film deliver the added hand protection needed for the most rigorous procedures
- **Confident grip:** These sterile gloves feature an external micro-textured finish for firmer, more confident instrument handling
- Effortless donning: Ansell's donning aid technology includes an internal polymer coating that facilitate easy dry or damp donning
- Chemo-tested: In addition, these long-cuff gloves are tested against ASTM D6978 and recommended for the safe handling of chemotherapy drugs\*

\*Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations.

#### Recommended For

- Damp donning
- Protection from latex Type I allergy and minimize chemical Type IV allergies in HCW's or patients
- · Tested for use with chemotherapy drugs

#### FEATURED TECHNOLOGIES









## GAMMEX<sup>®</sup> Innovation that Protects

Driven by innovation, GAMMEX<sup>®</sup> surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.



## GAMMEX® Non-Latex

SURGICAL GLOVES

Material	Neoprene (Polychloroprene)							
Powder Content	Powder-Free							
Color	Green							
Thickness	Standard							
Cuff Style	Straight with SUREFIT™ Technology							
External Glove Surface	Smooth with Micro-Textured Finish							
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology							
Allergy Prevention	Latex (Type I), Chemical (Type IV)							
Grip Level	Moderate							
Double Gloving Recommendation	Outer glove or underglove							
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)							
Freedom from Holes (Inspection level I)	0.65 AQL							
Sterile	Yes							
Product Standards	AS/NZS 4179, ASTM D3577, ASTM D6978, ASTM D7160, ASTM F1670, ASTM F739, EN 16523-1, EN 421:2010 EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282, JIS T9107							
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001, Korean GMP							
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case							
Shelf Life	3 years							
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.							
Vulcanization Chemical Accelerators	None							
ORDERING INFORMATION	N							
Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	340006055	340006060	340006065	340006070	340006075	340006080	340006085	34000609

#### Contact your Ansell representative for ordering or more information.

North America US Tel: 800 952 9916 CA Tel: 844 494 7854 Email: insidesalesus@ansell.com Sandel US Tel: 866 764 3327 Sandel email: sandel-cs@ansell.com

Central & South America Tel: +52(442) 296 20 50 Email: cslac@ansell.com Brazil Tel: +55-11-3356-3100 Email: latam.medical@ansell.com Europe, Middle East & Africa Tel: +32 (0) 2 528 74 00 Email: info@ansell.eu Asia Pacific Tel: +603 8310 6688 Email: apac.medical@ansell.com

India
Tet: +91 22 38124500
Email: info.india@ansell.com
Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com
China
Tet: +86(21) 3827 5005
Email: infochina@ansell.com