

# ENCORE<sup>®</sup> Latex Micro

S U R G I C A L   G L O V E S

## Latex, powder free micro surgical gloves offer exceptional tactile sensitivity and dexterity

- **Heightened dexterity:** Being 20% thinner than Ansell's standard powder-free latex gloves, ENCORE<sup>®</sup> Latex Micro offer exceptional level of tactile sensitivity and dexterity preferred for delicate procedures
- **Elevated ergonomic comfort:** Ansell's unique anatomic glove formers deliver surgical gloves designed for enhanced comfort and movement to help reduce hand fatigue
- **Double-gloving versatility:** This micro glove's brown color diminishes surgical glare while offering the wearer the option to use it as an inner glove when double gloving for easy breach detection
- **Enhanced donning aid:** A unique internal coating enables easy donning and doffing under dry or damp hand conditions, be it the first don or during intra-operative changes



CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

### Recommended For

- All surgeries, especially delicate procedures in ophthalmology, micro, and cardiovascular
- Damp donning
- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Prevention of eye fatigue

### FEATURED TECHNOLOGIES



## ENCORE<sup>®</sup> The Confidence to Perform

Dedicated to delivering optimized solutions engineered to provide a high degree of confidence to surgeons, O.R. nurses and surgical technicians seeking protection with a high level of dexterity, comfort and grip so they can perform at their best.

# ENCORE<sup>®</sup> Latex Micro

SURGICAL GLOVES

## PRODUCT INFORMATION

|   |   |
|---|---|
| Material                                | Natural Rubber Latex  |
| Powder Content                          | Powder-Free   |
| Color                                   | Brown   |
| Thickness                               | Micro   |
| Cuff Style                              | Beaded with SUREFIT™ Technology   |
| External Glove Surface                  | Smooth with Micro-Textured Finish   |
| Internal Glove Surface                  | Polymer Coated with DERMASHIELD™ Technology   |
| Allergy Prevention                      | None  |
| Grip Level                              | Moderate  |
| Double Gloving Recommendation           | Outer glove or underglove   |
| Tested For Use With Chemotherapy Drugs  | Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)  |
| Freedom from Holes (Inspection level I) | • In process control before packing: 0.65 (Inspection level GI) • Manufacturing final release: 0.65 (Inspection level GI)           |
| Sterile                                 | Yes   |
| Product Standards                       | ASTM D3577, ASTM D7160, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN 374-2 and -4, EN ISO 374-5:2016, ISO 10282 |
| Quality / Environmental Standards       | EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001  |
| Packaging                               | 50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case  |
| Shelf Life                              | 3 years   |
| Storage Instructions                    | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.                            |
| Vulcanization Chemical Accelerators     | Zinc Dibutyldithiocarbamate (ZDBC), Zinc Diethyldithiocarbamate (ZDEC)  |

## ORDERING INFORMATION

| Size         | 5.5       | 6         | 6.5       | 7         | 7.5       | 8         | 8.5       | 9         |
|--------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| Product Code | 330104055 | 330104060 | 330104065 | 330104070 | 330104075 | 330104080 | 330104085 | 330104090 |

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: insidesalesus@ansell.com  
 Sandel US Tel: 866 764 3327  
 Sandel email: sandel-cs@ansell.com

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: cslac@ansell.com  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: [latam.medical@ansell.com](mailto:latam.medical@ansell.com)

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: info@ansell.eu  
**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**India**  
 Tel: +91 22 38124500  
 Email: info.india@ansell.com  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com  
**China**  
 Tel: +86(21) 3827 5005  
 Email: infochina@ansell.com