

# GAMMEX® Latex Moisturizing

S U R G I C A L   G L O V E S

## Latex powder-free surgical gloves, featuring a built-in moisturizer to rehydrate the skin

- **Soft, supple skin:** GAMMEX® Latex Moisturizing surgical gloves feature a proprietary skin moisturizing technology, enabling a water-based coating with glycerol and dimethicone to retain moisture and rehydrate dry skin, improving the wearer's skin condition
- **Exceptional comfort:** These gloves are also designed using anatomic formers for added comfort and fit
- **Enhanced grip:** With a smooth finish and micro-textured external surface around the palm and fingers, these surgical gloves deliver excellent grip under dry or damp conditions
- **Easy donning:** Ansell's donning aid technology facilitates an easy first don, whether the wearer's hands are dry or damp
- **Optimal allergy prevention:** Additional cleaning processes ensure low residues of latex proteins and chemical accelerators, minimizing the risk of developing latex Type I and chemical Type IV allergies and sensitivities



**CAUTION:** Safe use of these gloves by latex-sensitized individuals has not been established – this product contains natural rubber latex, which may cause allergic reactions

### Recommended For

- Damp donning
- Extra comfort during long, sustained procedures
- Use as an underglove when double gloving

### FEATURED TECHNOLOGIES



### GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX<sup>®</sup> Latex Moisturizing

## SURGICAL GLOVES

### PRODUCT INFORMATION

|                                         |                                                                                                                    |
|-----------------------------------------|--------------------------------------------------------------------------------------------------------------------|
| Material                                | Natural Rubber Latex                                                                                               |
| Powder Content                          | Powder-Free                                                                                                        |
| Color                                   | Natural                                                                                                            |
| Thickness                               | Standard                                                                                                           |
| Cuff Style                              | Straight with SUREFIT™ Technology                                                                                  |
| External Glove Surface                  | Smooth with Micro-Textured Finish                                                                                  |
| Internal Glove Surface                  | Coated with HYDRASOFT™ Technology                                                                                  |
| Allergy Prevention                      | None                                                                                                               |
| Grip Level                              | Moderate                                                                                                           |
| Double Gloving Recommendation           | Underglove                                                                                                         |
| Tested For Use With Chemotherapy Drugs  | No, not tested in accordance with ASTM D6978                                                                       |
| Freedom from Holes (Inspection level I) | 0.65 AQL                                                                                                           |
| Sterile                                 | Yes                                                                                                                |
| Product Standards                       | AS/NZS 4179, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282, JIS T9107 |
| Quality / Environmental Standards       | EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001                                                           |
| Packaging                               | 50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case                                               |
| Shelf Life                              | 3 years                                                                                                            |
| Storage Instructions                    | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.           |
| Vulcanization Chemical Accelerators     | Xanthogen Polysulphides, Sodium di-n-butyldithiocarbamate (SDBC), Zinc Diethyldithiocarbamate (ZDEC)               |

### ORDERING INFORMATION

| Size         | 5.5       | 6         | 6.5       | 7         | 7.5       | 8         | 8.5       | 9         |
|--------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| Product Code | 330049055 | 330049060 | 330049065 | 330049070 | 330049075 | 330049080 | 330049085 | 330049090 |

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: [insidesalesus@ansell.com](mailto:insidesalesus@ansell.com)  
 Sandel US Tel: 866 764 3327  
 Sandel email: [sandel-cs@ansell.com](mailto:sandel-cs@ansell.com)

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: [csiac@ansell.com](mailto:csiac@ansell.com)  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: [latam.medical@ansell.com](mailto:latam.medical@ansell.com)

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: [info@ansell.eu](mailto:info@ansell.eu)  
**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: [apac.medical@ansell.com](mailto:apac.medical@ansell.com)

**India**  
 Tel: +91 22 38124500  
 Email: [info.india@ansell.com](mailto:info.india@ansell.com)  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: [protection@ap.ansell.com](mailto:protection@ap.ansell.com)  
**China**  
 Tel: +86(21) 3827 5005  
 Email: [infochina@ansell.com](mailto:infochina@ansell.com)