

# GAMMEX<sup>®</sup> Latex EnLite

S U R G I C A L   G L O V E S

## Latex powder-free surgical gloves, for reliable hand protection and a comfortable fit

- **Elevated comfort:** GAMMEX<sup>®</sup> Latex EnLite™ latex surgical gloves are made using unique anatomic formers, ensuring an excellent fit and superb comfort
- **Effortless donning:** Thanks to a chlorination treatment applied to its internal surface, these latex gloves are easy to don
- **Secure hold:** GAMMEX<sup>®</sup> Latex EnLite™ sterile gloves feature straight cuffs, ensuring a secure hold to the gown
- **Optimal allergy prevention:** An exclusive manufacturing process ensures very low levels of latex proteins and chemical accelerator residues, minimizing the risk of developing latex Type I and chemical Type IV allergies and sensitivities
- **Secure grip:** These gloves' smooth surface, with a micro-textured finish, enable confident instrument handling
- **Chemo-drugs tested:** Tested for safe use of chemo-drugs in accordance with ASTM D6879 requirements\*



**CAUTION:** Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

\*Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

### Recommended For

- Brief, low-risk examination procedures
- Extra comfort during long, sustained procedures
- Sterile procedures in the O.R.
- Use as an outer glove when double gloving

### FEATURED TECHNOLOGIES



## GAMMEX<sup>®</sup> Innovation that Protects

Driven by innovation, GAMMEX<sup>®</sup> surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX<sup>®</sup> Latex EnLite

SURGICAL GLOVES

## PRODUCT INFORMATION

<b>Material</b>	Natural Rubber Latex
<b>Powder Content</b>	Powder-Free
<b>Color</b>	Natural
<b>Thickness</b>	Standard
<b>Cuff Style</b>	Straight
<b>External Glove Surface</b>	Smooth with Micro-Textured Finish
<b>Internal Glove Surface</b>	Chlorinated
<b>Allergy Prevention</b>	None
<b>Grip Level</b>	Moderate
<b>Double Gloving Recommendation</b>	Outer glove
<b>Tested For Use With Chemotherapy Drugs</b>	No, not tested in accordance with ASTM D6978
<b>Freedom from Holes (Inspection level I)</b>	0.65 AQL
<b>Sterile</b>	Yes
<b>Product Standards</b>	EN 420:2003 + A1:2009, EN 455 1-4, ISO 10282, JIS T9107
<b>Quality / Environmental Standards</b>	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
<b>Packaging</b>	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
<b>Shelf Life</b>	3 years
<b>Storage Instructions</b>	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
<b>Vulcanization Chemical Accelerators</b>	Zinc Diethyldithiocarbamate (ZDEC)

## ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
<b>Product Code</b>	330061055	330061060	330061065	330061070	330061075	330061080	330061085	330061090

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: insidesalesus@ansell.com  
 Sandel US Tel: 866 764 3327  
 Sandel email: sandel-cs@ansell.com

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: cslac@ansell.com  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: [latam.medical@ansell.com](mailto:latam.medical@ansell.com)

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: info@ansell.eu  
**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**India**  
 Tel: +91 22 38124500  
 Email: info.india@ansell.com  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com  
**China**  
 Tel: +86(21) 3827 5005  
 Email: infochina@ansell.com