

# GAMMEX® Non-Latex PI Underglove

S U R G I C A L   G L O V E S

## Soft, skin-friendly polyisoprene surgical gloves, for heightened tactile sensitivity

- **Superior comfort and protection:** Ansell's proprietary soft yet strong formulation delivers heightened comfort and hand protection when double gloving
- **Easy breach detection:** Green in color, these sterile gloves enable fast and easy breach detection when worn with a cream or white outer glove
- **Optimal allergy protection:** Being non-latex, these polyisoprene surgical gloves eliminate the risk of developing a latex allergy while their skin-friendly formulation helps minimize sensitization to developing a Type IV chemical allergy as chemical accelerators known to cause skin allergies have been eliminated
- **Effortless donning and doffing:** A unique internal coating enables easy donning and doffing under dry or damp hand conditions



## Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Easy identification of a glove breach
- Protection from chemical accelerators such as thiurams, carbamates and mercaptos known to cause chemical Type IV allergies
- Having passed the Mouse Embryo Assay (MEA) test, may be used for In Vitro Fertilization (IVF) procedures

## FEATURED TECHNOLOGIES



## GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.



# GAMMEX® Non-Latex PI Underglove

## SURGICAL GLOVES

### PRODUCT INFORMATION

|                                         |                                                                                                                                             |
|-----------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| Material                                | PI (Polyisoprene)                                                                                                                           |
| Powder Content                          | Powder-Free                                                                                                                                 |
| Color                                   | Green                                                                                                                                       |
| Thickness                               | Micro                                                                                                                                       |
| Cuff Style                              | Beaded with SUREFIT™ Technology                                                                                                             |
| External Glove Surface                  | Smooth with Micro-Textured Finish                                                                                                           |
| Internal Glove Surface                  | Polymer Coated with DERMASHIELD™ Technology                                                                                                 |
| Allergy Prevention                      | Latex (Type I)                                                                                                                              |
| Grip Level                              | Low                                                                                                                                         |
| Double Gloving Recommendation           | Underglove                                                                                                                                  |
| Tested For Use With Chemotherapy Drugs  | Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)                                                                          |
| Freedom from Holes (Inspection level I) | 0.65 AQL                                                                                                                                    |
| Sterile                                 | Yes                                                                                                                                         |
| Product Standards                       | AS/NZS 4179, ASTM D3577, ASTM D6978, ASTM D7160, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282 |
| Quality / Environmental Standards       | EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001                                                                                    |
| Packaging                               | 50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case                                                                        |
| Shelf Life                              | 3 years                                                                                                                                     |
| Storage Instructions                    | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.                                    |
| Vulcanization Chemical Accelerators     | Xanthates, Zinc Dialkyldithiocarbamate                                                                                                      |

### ORDERING INFORMATION

| Size         | 5.5       | 6         | 6.5       | 7         | 7.5       | 8         | 8.5       | 9         |
|--------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| Product Code | 340058055 | 340058060 | 340058065 | 340058070 | 340058075 | 340058080 | 340058085 | 340058090 |

Contact your Ansell representative for ordering or more information.

**North America**  
US Tel: 800 952 9916  
CA Tel: 844 494 7854  
Email: [insidesalesus@ansell.com](mailto:insidesalesus@ansell.com)  
Sandel US Tel: 866 764 3327  
Sandel email: [sandel-cs@ansell.com](mailto:sandel-cs@ansell.com)

**Central & South America**  
Tel: +52(442) 296 20 50  
Email: [cslac@ansell.com](mailto:cslac@ansell.com)  
**Brazil**  
Tel: +55-11-3356-3100  
Email: [latam.medical@ansell.com](mailto:latam.medical@ansell.com)

**Europe, Middle East & Africa**  
Tel: +32 (0) 2 528 74 00  
Email: [info@ansell.eu](mailto:info@ansell.eu)  
**Asia Pacific**  
Tel: +603 8310 6688  
Email: [apac.medical@ansell.com](mailto:apac.medical@ansell.com)

**India**  
Tel: +91 22 38124500  
Email: [info.india@ansell.com](mailto:info.india@ansell.com)  
**Australia & New Zealand**  
Tel: +61 3 9270 7270  
Email: [protection@ap.ansell.com](mailto:protection@ap.ansell.com)  
**China**  
Tel: +86(21) 3827 5005  
Email: [infochina@ansell.com](mailto:infochina@ansell.com)

**ansell.com**

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking).  
© 2025 Ansell Limited. All Rights Reserved.