

GAMMEX® Non-Latex PI Textured

SURGICAL GLOVES

Skin-friendly polyisoprene surgical gloves featuring a textured external finish for superb grip

- **Enhanced grip:** These synthetic surgical gloves' textured finish enhances instrument grip, for secure handling in wet environments, especially during high-risk procedures
- **Easy double gloving:** Best used as an outer glove when double gloving, Ansell's proprietary donning aid technology includes an inner polymer coating to facilitate donning and doffing during intra-operative changes
- **Optimized allergenic prevention:** These sterile polyisoprene gloves eliminate the risk of developing a latex allergy and help prevent sensitization to allergic contact dermatitis, as they are not made with chemical accelerators known to be allergens
- **Chemo-tested:** GAMMEX® Non-Latex PI Textured gloves are tested for handling with a range of cytotoxic drugs per ASTM D6978**



**Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Enhanced grip in wet environments
- Protection from chemical accelerators such as thiurams, carbamates and mercaptos known to cause chemical Type IV allergies
- Use as an outer glove when double gloving
- Having passed the Mouse Embryo Assay (MEA) test, may be used for In Vitro Fertilization (IVF) procedures



FEATURED TECHNOLOGIES



GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.



GAMMEX[®] Non-Latex PI Textured

SURGICAL GLOVES

PRODUCT INFORMATION

Material	PI (Polyisoprene)
Powder Content	Powder-Free
Color	White
Thickness	Standard
Cuff Style	Beaded with SUREFIT [™] Technology
External Glove Surface	Textured Fingers and Palm
Internal Glove Surface	Polymer Coated with DERMASHIELD [™] Technology
Allergy Prevention	Latex (Type I)
Grip Level	High
Double Gloving Recommendation	Outer glove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Freedom from Holes (Inspection level I)	0.4 AQL
Sterile	Yes
Product Standards	AS/NZS 4179, ASTM D3577, ASTM D6978, ASTM D7160, ASTM F739, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282, JIS T9107
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001, Korean GMP
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Xanthates, Zinc Dialkyldithiocarbamate

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	340054055	340054060	340054065	340054070	340054075	340054080	340054085	340054090

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 844 494 7854
Email: insidesalesus@ansell.com
Sandel US Tel: 866 764 3327
Sandel email: sandel-cs@ansell.com

Central & South America
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil
Tel: +55-11-3356-3100
Email: latam.medical@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu
Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

India
Tel: +91 22 38124500
Email: info.india@ansell.com
Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com
China
Tel: +86(21) 3827 5005
Email: infochina@ansell.com

ansell.com

Ansell, [®] and [™] are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking.
© 2025 Ansell Limited. All Rights Reserved.