

# GAMMEX<sup>®</sup> Latex Underglove

S U R G I C A L   G L O V E S

## Green latex under gloves for easy breach detection and moisturizing hand protection when double gloving

- **Optimized skin hydration:** These powder-free latex under gloves feature an odorless, water-based skin moisturizer containing glycerol and dimethicone\* to help retain moisture and rehydrate skin throughout the duration of use
- **Added safety:** When contrasted with a cream or white-colored outer glove, these green under gloves offer quick and easy breach detection for glove wearer's added protection
- **Effortless donning and doffing:** The internal polymer coating of the GAMMEX<sup>®</sup> Latex Underglove helps to facilitate effortless dry and damp hand donning, whether it be the first don or during intra-operative changes
- **Elevated comfort:** The GAMMEX<sup>®</sup> Latex Underglove is sterile and made using proprietary anatomic-shaped formers and a soft formulation for ergonomic fit and comfort to help minimize hand fatigue during rigorous procedures

\*built-in moisturizing agents

**CAUTION:** Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

### Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Damp donning
- Easy identification of a glove breach
- Use as an underglove when double gloving

### FEATURED TECHNOLOGIES



### GAMMEX<sup>®</sup> Innovation that Protects

Driven by innovation, GAMMEX<sup>®</sup> surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX<sup>®</sup> Latex Underglove

SURGICAL GLOVES

## PRODUCT INFORMATION

<b>Material</b>	Natural Rubber Latex
<b>Powder Content</b>	Powder-Free
<b>Color</b>	Green
<b>Thickness</b>	Standard
<b>Cuff Style</b>	Straight with SUREFIT™ Technology
<b>External Glove Surface</b>	Smooth with Micro-Textured Finish
<b>Internal Glove Surface</b>	Coated with HYDRASOFT™ Technology
<b>Allergy Prevention</b>	None
<b>Grip Level</b>	Low
<b>Double Gloving Recommendation</b>	Underglove
<b>Tested For Use With Chemotherapy Drugs</b>	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
<b>Freedom from Holes (Inspection level I)</b>	0.65 AQL
<b>Sterile</b>	Yes
<b>Product Standards</b>	AS/NZS 4179, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282
<b>Quality / Environmental Standards</b>	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001, Korean GMP
<b>Packaging</b>	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
<b>Shelf Life</b>	3 years
<b>Storage Instructions</b>	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
<b>Vulcanization Chemical Accelerators</b>	Xanthogen Polysulphides, Sodium di-n-butylthiocarbamate (SDBC), Zinc Diethyldithiocarbamate (ZDEC)

## ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
<b>Product Code</b>	330050055	330050060	330050065	330050070	330050075	330050080	330050085	330050090

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: insidesalesus@ansell.com  
 Sandel US Tel: 866 764 3327  
 Sandel email: sandel-cs@ansell.com

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: cslac@ansell.com  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: [latam.medical@ansell.com](mailto:latam.medical@ansell.com)

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: info@ansell.eu  
**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**India**  
 Tel: +91 22 38124500  
 Email: info.india@ansell.com  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com  
**China**  
 Tel: +86(21) 3827 5005  
 Email: infochina@ansell.com