

MEDI-GRIP® Latex Powder-Free

S U R G I C A L G L O V E S

Latex, powder-free surgical gloves with a micro-textured finish for secure instrument handling

- **Enhanced grip:** The unique external surface treatment and micro-textured finish of MEDI-GRIP® Latex Powder-Free surgical gloves enables the glove wearer to securely grip medical instruments for safe handling under dry or damp conditions
- **Elevated comfort:** The unique anatomic formers and soft latex formulation of these powder-free sterile gloves ensure excellent fit and feel, offering prolonged comfort during long procedures and low risk examinations
- **Excellent dexterity:** Their soft, malleable formulation ensures dexterity to support easy manipulation of surgical instruments being used
- **Chemo-tested:** Tested to ASTM D6978, MEDI-GRIP® Latex Powder-Free gloves are safe to use with a wide range of chemotherapy drugs*



CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

*Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

Recommended For

- Sterile procedures in the O.R.
- Extra comfort during long, sustained procedures
- Brief, low-risk examination procedures
- Use as an outer glove when double gloving

FEATURED TECHNOLOGIES



MEDI-GRIP®
Everyday Protection you can
Depend on

Trusted around the world to deliver the every day protection surgeons, O.R. nurses and surgical technicians can depend on to perform their daily tasks safely and comfortably.

MEDI-GRIP[®] Latex Powder-Free

SURGICAL GLOVES

PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder-Free
Color	Natural
Thickness	Standard
Cuff Style	Beaded
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated
Allergy Prevention	None
Grip Level	Moderate
Double Gloving Recommendation	Outer glove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	100 pairs per box; 3 boxes per carton/case; 300 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC)

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	354101	354102	354103	354104	354105	354106	354107	354108

Contact your Ansell representative for ordering or more information.

North America
 US Tel: 800 952 9916
 CA Tel: 844 494 7854
 Email: insidesalesus@ansell.com
 Sandel US Tel: 866 764 3327
 Sandel email: sandel-cs@ansell.com

Central & South America
 Tel: +52(442) 296 20 50
 Email: cslac@ansell.com
Brazil
 Tel: +55-11-3356-3100
 Email: latam.medical@ansell.com

Europe, Middle East & Africa
 Tel: +32 (0) 2 528 74 00
 Email: info@ansell.eu
Asia Pacific
 Tel: +603 8310 6688
 Email: apac.medical@ansell.com

India
 Tel: +91 22 38124500
 Email: info.india@ansell.com
Australia & New Zealand
 Tel: +61 3 9270 7270
 Email: protection@ap.ansell.com
China
 Tel: +86(21) 3827 5005
 Email: infochina@ansell.com