

MICRO-TOUCH® Nitrile Accelerator-Free

H i g h T o u c h

Nitrile, powder-free medical exam glove with advanced allergy prevention


- **Comfortable fit:** This non-sterile, accelerator-free nitrile glove offers a comfortable fit and feel
- **Barrier protection:** This glove provides a high level of barrier protection against punctures and excellent resistance from a broad range of commonly used chemicals
- **Easy donning:** Chlorinated inner surface allows for easy donning
- **Reduced allergy risk:** Nitrile material prevents latex Type I allergies — 100% thiuram-free formulation significantly reduces the risk of Type IV allergies
- **Secure grip:** Micro-textured fingers ensure a secure and safe grip in wet or dry conditions
- **Specialized protection:** Tested and approved for use with chemotherapy drugs*

*see product packaging for details

Recommended For

- Tested for use with chemotherapy drugs
- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Protection from latex Type I allergy and minimize chemical Type IV allergies in HCW's or patients
- Standard patient examinations

FEATURED TECHNOLOGIES



MICRO-TOUCH® Everyday, Everywhere Protection

MICRO-TOUCH® premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH® Nitrile Accelerator-Free

High Touch

PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Blue	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.08 / 3.1
Finger Thickness (mm/mil)	0.12 / 4.7	Allergy Prevention	Latex (Type I), Chemical (Type IV)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	240 / 9.5	Product Segmentation	High Touch

ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	700101	700102	700103	700104	700105

PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

CE 0493, Category III, EN 16523-1, EN 421:2010, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 9001, ISO 10993-10:2002, ISO 13485

Vulcanization Chemical Accelerators

Accelerator Free

Only a very small number of users may be sensitive to this ingredient(s) and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email:
insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email:
luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email:
apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email:
protection@ap.ansell.com