

Sterile Cleanroom Pen

BioClean Permaflow™ Cleanroom Pens are sterile with fade and water resistant super permanent ink. With alcohol resistant labeling for easy cleaning within a critical environment, they are ISO Class 4 & EU GMP Grade A compatible. Available in a range of colours.

KEY FEATURES & BENEFITS

- Autoseal ink prevents drying
- Alcohol resistant labelling
- Processed to ensure ISO Class 4 & EU GMP Grade A compatibility
- Super permanent ink Black, Blue, Red or Green
- Fade and water resistant
- Quick drying
- Non-toxic

Industries

- Controlled and Critical Environments
- Laboratory and Research
- Production and Manufacturing



S-BPFP

Sterile Cleanroom Pen

TECHNICAL DATA SHEET

	Product Information
Material	Plastic
Color	Blue, Black
Manufacturing/QMS Audit Standards	ISO 9001:2015
Regulatory/Standard Compliance	ISO 9001:2015
Packaging	One pen per sealed inner PE bag; 10 inner bags per sealed outer PE bag; 20 outer bags per lined carton (200 pens)
Storage	Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.
Sterilization Method	GAMMA irradiation (25 kGy)
Sterility Assurance Level	10 ⁻⁶
Cleanroom Class	Class 10/ISO 4
Shelf Life	Five (5) years from date of manufacture.



S-BPFP

Sterile Cleanroom Pen

PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE

For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
US T: +1 800 800 0444
US F: +1 800 800 0445
CA T: +1-800-363-8340

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

UK

Ansell Nitritex
T: +44 1638 663338
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending; www.ansell.com/patentmarking © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

