

# MICRO-TOUCH<sup>®</sup> DermaClean<sup>®</sup>

Essential Series

## Latex, powder-free medical exam glove with excellent grip properties

- **Secure grip:** Thanks to its textured outer surface, this non-sterile glove has excellent grip levels in wet and dry conditions, and is ideal for healthcare providers who handle instrumentation, test tubes or infusion equipment
- **Comfortable fit:** Ansell's natural rubber latex formulation ensures flexibility in a comfortable, well-fitting glove
- **Easy donning:** The smooth, chlorinated inner surface facilitates ease of donning

CAUTION: This product contains natural rubber latex, which may cause allergic reactions

### Recommended For

- Brief, low-risk examination procedures



## MICRO-TOUCH<sup>®</sup> Everyday, Everywhere Protection

MICRO-TOUCH<sup>®</sup> premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH® DermaClean®

Essential Series

## PRODUCT INFORMATION

|   |                      |                                    |  |
|---|----------------------|------------------------------------|--|
| Material                                | Natural Rubber Latex | Not Made From Natural Rubber Latex | No   |
| Color                                   | White                | Cuff Length                        | Standard   |
| Powder Content                          | Powder-Free          | External Glove Surface             | Fully Textured   |
| Freedom from Holes (Inspection level I) | 1.5 AQL              | Finger Thickness (mm/mil)          | 0.16 / 6.3   |
| Palm Thickness (mm/mil)                 | 0.14 / 5.5           | Available Sizes                    | XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10) |
| Tested For Use With Chemotherapy Drugs  | No                   | Silicone Free                      | Yes  |
| Sterile                                 | No                   | Antistatic                         | Not Tested   |
| Glove Length (mm/inches)                | 240 / 9.5            | Product Segmentation               | Essential Series   |

## ORDERING INFORMATION

| Size         | XS (5-5.5) | S (6-6.5) | M (7.5 - 8) | L (8.5 - 9) | XL (9.5 - 10) |
|--------------|------------|-----------|-------------|-------------|---------------|
| Product Code | 3180-1     | 3181-1    | 3182-1      | 3183-1      | 3184-1        |

## PACKAGING AND STORAGE

|                      |  |
|----------------------|--|
| Packaging            | 100 gloves per dispenser 10 dispensers per case 1000 gloves per case *80 gloves per dispenser (Size XL) *800 gloves per case (Size XL) |
| Shelf Life           | 3 Years  |
| Storage Instructions | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.                               |

## STANDARDS AND CERTIFICATIONS

CA Prop 65, CE Number, Medical Device Regulation (EU) 2017/745 Class I, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 13485, ISO 10993-5, ISO 9001, ISO 11193-1:2008, ISO 10993-1, ISO 16604, ISO 10993-10, ISO 10993-10:2002, Food Type: Rubber

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email: insidesalesus@ansell.com

Latin America & Caribbean  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

Asia Pacific  
Tel: +603 8310 6688  
Email: apac.medical@ansell.com

Australia & New Zealand  
Tel: +61 3 9270 7270  
Email: protection@ap.ansell.com