

# MICRO-TOUCH® Nitrile Accelerator-Free

High Touch

## Nitrile, powder-free medical exam glove with advanced allergy prevention

- **Comfortable fit:** This non-sterile, accelerator-free nitrile glove offers a comfortable fit and feel
- **Barrier protection:** This glove provides a high level of barrier protection against punctures and excellent resistance from a broad range of commonly used chemicals
- **Easy donning:** Chlorinated inner surface allows for easy donning
- **Reduced allergy risk:** Manufactured without accelerators, helping to protect users from Type IV chemical allergies
- **Secure grip:** Micro-textured fingers ensure a secure and safe grip in wet or dry conditions
- **Specialized protection:** Tested and approved for use with chemotherapy drugs\*




\*see product packaging for details

### Recommended For

- Tested for use with chemotherapy drugs
- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Protection from latex Type I allergy and minimize chemical Type IV allergies in HCW's or patients
- Standard patient examinations

### FEATURED TECHNOLOGIES



### MICRO-TOUCH® Everyday, Everywhere Protection

MICRO-TOUCH® premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH® Nitrile Accelerator-Free

High Touch

## PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Blue	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Finger Thickness (mm/mil)	0.12 / 4.7
Palm Thickness (mm/mil)	0.08 / 3.1	Allergy Prevention	Latex (Type I), Chemical (Type IV)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	240 / 9.5	Product Segmentation	High Touch

## ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	3130230XS	31302300S	31302300M	31302300L	3130230XL

## PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

## STANDARDS AND CERTIFICATIONS

CE 0493, Category III, Medical Device Regulation (EU) 2017/745 Class I, EN 16523-1, EN 421:2010, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10993-10:2002, ISO 9001, ISO 13485, UKCA

Contact your Ansell representative for ordering or more information.

**North America**  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email: insidesalesus@ansell.com

**Latin America & Caribbean**  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email: luvas.medicas@ansell.com

**Europe, Middle East & Africa**  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

**Asia Pacific**  
Tel: +603 8310 6688  
Email: apac.medical@ansell.com

**Australia & New Zealand**  
Tel: +61 3 9270 7270  
Email: protection@ap.ansell.com